
(For official use only)

(Govt. of India)
(Ministry of Railways)

HAND BOOK ON TESTING SCHEME FOR BIO-TOILET EFFLUENT & BACTERIAL CULTURE (INOCULUM)

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अभ्यास RDSO
रेल अग्रदूत Transforming Railways

Indian Railways
Centre for Advanced Maintenance Technology

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Testing Scheme for Bio-Toilet Effluent and Bacteria culture (Inoculum) on Indian Railways

On the basis of MOU between Indian Railways and DRDO, the bio toilet trial on IR is going on. For intensive monitoring of the Bio-Toilet system, a test scheme has been developed for implementation on IR after formal approval of RDSO/LKO. As per Railway Boards guidelines, the bio-toilet tanks i.e. Bio digesters are being fitted in the new coaches being turned out from RCF. At present, the initial charging of bacterial culture in the coaches is being done by DRDE/ Gwalior by providing Inoculum in the barrels. The quality check of the product and bio toilet effluent is also done according to ALPHA testing methods in DRDE by collecting samples from coaching depot / Gwalior,

As per Railway Board's guidelines, a 100 M3 capacity Inoculum generation plant has been developed in the Motibagh workshop in S.E.C. Railway, Nagpur and in future, similar plants are proposed to be developed in ICF and RCF also. Hence, it was felt necessary to develop lab sample testing facilities in Railways for doing minor testings in Railway premises/Railway Labs/ production units & major testings in Govt. approved labs. To cater the requirement of doing these tests efficiently, it has also been proposed to provide suitable training to Railway staff in liaison with DRDE so that proper monitoring can be ensured at appropriate level in the zonal Railways. Some of the minor tests of effluent of bio-toilets have been proposed to be carried out by coaching depot staff for which procedure has been explained in detail in this booklet.

The testing scheme for effluent & Bacterial culture(Inoculum) has been tentatively finalized by DRDE based on the past experience of 12 months trial, which should be carried out in the existing Railway Labs/ diesel sheds/ coaching depots. The necessary record of monitoring and testing should also be maintained & kept in the concerned depots. The critical biological tests which can not be carried out in the railway labs, such samples may be got tested in govt. labs/ DRDE. The critical tests which should be carried out in Govt. labs only, has been marked clearly against the test procedure in the booklet for ready reference

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Government of India-Ministry of Railways
Research Designs & Standards Organisation
Lucknow - 226 011
DID (0522) 2450115
DID (0522) 2465310

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Sub: Issue of trial/test scheme for field trial of bio toilets (DRDE Technology) on Indian Railway Passenger Coaches.

Ref: 1. Item No. 12 of the MOM of 7th JWG held at RDSO on 01.10.2012.
2. This office letter of even No. dated 21.06.2012

Please find enclosed revised trial/test scheme number RDSO/2010/CG/TS10 Rev-04 incorporating the sampling procedure for the effluent of IR DRDO Bio Toilets along with the list of Government accredited test labs, based on the inputs given by DRDE Gwalior for field trial of IR-DRDE type Bio-toilets. It is for your kind information and necessary action.

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Government of India
Ministry of Defence
Defence R&D Organization
Defence R&D Establishment
Jhansi Road, Gwalior- 474 002
(MP) INDIA

Gram : DEFRES
Telex : 0786- 212
Phone : 0751-2340354, 2340245,
2341848 Ext 270
Fax : 91-751-2341148
Email : drde@sancharnet.in
Date : 26.03.2012

Ref BT/ BKB/ 002/biodig/2011-12

To,

The Executive Director,
Centre for Advanced maintenance technology (CAMTECH)
Maharajpur, Gwalior – 474 005

Sub: IR- DRDO Bio- toilets.

Ref : IRCAMTECH/GWL/M/Bio – Toilets dated 27.02.2012

Please find enclosed herewith the test procedures developed for bacterial inoculum facility and effluent discharge parameters.

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Testing Scheme for Bio-Toilet Effluent and Bacterial culture (Inoculum)

For proper monitoring of Bio-Toilet system, timely testing of samples should be ensured at appropriate level in the zonal Railways. Following 04 lab tests for effluent of bio-toilets have been proposed to be carried out by coaching depot staff for which procedure has been explained in detail in this hand book.

The Testing Scheme for Effluent & Bacterial culture(Inoculum) has been tentatively finalized by DRDE/GWL based on the past experience of 12 months trial, which should be carried out in the existing Railway Labs/ diesel sheds/ coaching depots. The necessary record of monitoring and testing should also be maintained & kept in the Coaching depots. The critical biological tests like COD test and Fecal Coli Forms count & MPN count for methanogens which can not be carried out in the railway labs, should be got tested in govt. labs/ DRDE. For more details & guidelines on lab sample testing, refer RDSO Doc No. RDSO/2010/CG/TS-10 Revision 04 dated 19.11.2012.

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(A) Testing Scheme for Bio-Toilet Effluent in Coaching Depots.

Following tests and parameters have been recommended for effluent discharged after bio degradation from the bio toilets of coaches in Indian Railways.

1. pH Value Test:

To measure pH value of the effluent of bio toilets to ensure environmental safety.


Table top pH meter

Portable pH meter

pH indicator strips

S. No.	Description	Details
1.	Purpose of test	To measure pH value of the effluent of bio toilets to ensure environmental safety.
2.	Target value	6 – 9 pH
3.	Equipments required	Table top pH meter / Portable pH meter / pH indicator strips & magnetic stirrer
4.	Consumables	pH calibration buffers (4.0, 7.0, 10.0), and magnetic stirrer bars
5.	Quantity of sample	50-100 ml
6.	Electricity requirement	Yes
7.	Frequency of sampling	90 days
8.	Testing spot	Railway laboratory
9.	Staff required	1

Procedure:

- First of all, take 02 ltrs. Effluent sample in the bottle from bio-toilet tank.
- Take 50 – 100 ml of mixed effluent sample in a beaker.
- Put a magnetic bar and keep the beaker on a magnetic stirrer and switch “on” the magnetic stirrer to mix it continuously.
- Wash the electrode and Temperature compensation rod with distilled water and wipe it with tissue paper.
- Put the electrode and automatic temperature compensation rod into the sample and keep it until stable reading appears in the displays; note the reading.
- Discard the sample and wash the electrode and automatic temperature compensation rod with distilled water and wipe it with tissue paper.
- Keep the electrode back in the container.

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pH Value Test (Simplified Procedure):

To measure pH value of the effluent of bio toilets to ensure environmental safety.

S. No.	Description	Details
01	Purpose of test	To measure pH value of the effluent of bio toilets to ensure environmental safety.
02	Target value	6 – 9 pH
03	Equipments required	Portable pH meter
04	Other equipments	pH indicator paper of all ranges
05	Consumables	pH calibration buffers (4.0, 7.0, 10.0)
06	Quantity of sample	50-100 ml
07	Electricity requirement	Yes
08	Frequency of sampling	90 days
09	Testing spot	Railway laboratory
10	Staff required	01

Procedure:

- First of all, take 02 ltrs. Effluent sample in the bottle from bio-toilet tank.
- Take 50 – 100 ml of mixed effluent sample in a beaker.
- Put pH meter on the stand .
- Mix the liquid properly.
- Put the electrode of portable pH meter in the beaker until stable reading appears in the displays;
- Note the reading from display of Ph meter.
- Discard the sample and wash the electrode of pH meter with distilled water and wipe it with tissue paper.
- Keep the electrode back in the container.

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2. Total Solids (TS) Test:

This test should be carried out to estimate amount of total solids available in the effluent.

Hot Air Oven

Weighing balance

Silica crucible

S. N	Description	Details
1.	Purpose of test	To estimate amount of total solids in the effluent.
2.	Set Target	<750 mg /100 ml
3.	Equipments required	Electronic weighing balance, pipettes, Silica crucible, Hot air oven, desiccators.
4.	Consumables	Self indicating silica gel.
5.	Quantity of sample	25 ml
6.	Electricity requirement	Yes
7.	Frequency of sampling	90 days
8.	Testing spot	Railway laboratory
9.	Staff required	1

Procedure:

- Heat an empty crucible and clean silica at 103 – 105°C for 1 hour in a hot air oven.
- Cool in desiccator to room temperature and take the initial weight.
- Pipette a measured volume of well mixed sample (25 ml).
- Keep the silica crucible in a hot air oven at 103 – 105°C for 1 hour; keep till the water gets dried or constant weight is achieved.
- Remove the silica crucible and keep it in desiccator until it reaches room temperature.
Note the final weight

Calculations:

$$\text{mg total solids/ 100 mL} = (A - B) \times 100 \times 1000 / \text{Volume of sample (mL)}$$

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Silica crucible

Filter flask

S. N	Description	Details
1.	Purpose of test	To estimate amount of total dissolved solids in the effluent.
2.	Target value	< 350 mg /100 ml
3.	Equipments required	Electronic weighing balance, pipettes, Silica crucibles, Hot air oven, desiccators, filter assembly, vacuum pump
4.	Consumable	Self indicating silica gel. Whatman Glass wool filters
5.	Quantity of sample	25 ml
6.	Electricity requirement	Yes
7.	Frequency of sampling	90 days
8.	Testing spot	Railway laboratory
9.	Staff required	1

Procedure:

1. Heat an empty and clean silica crucible at $180 \pm 2^{\circ}\text{C}$ for 1 to /2 hour in a hot air oven, cool it in desiccator to room temperature and weight.
2. Insert disc and filter assembly, apply vacuum and wash disc with three successive 20-mL volumes of reagent grade water.
3. Dry it in a hot air oven at $180 \pm 2^{\circ}\text{C}$ for 1 hour in an oven. Store in desiccators until needed. Weigh immediately before use.
4. Stir sample with a magnetic stirrer and pipette a measured volume (25 ml) into a glass fiber filter with applied vacuum.
5. Wash with three successive 10-mL volume of distilled water, allowing complete drainage between washings.
6. Transfer total filtrate (with washings) to a pre-weighed clean silica crucible.
7. Dry it in a hot air oven at $180 \pm 2^{\circ}\text{C}$ for 1 hour in an oven. Cool it in desiccators to room temperature. Note the final weight

Calculations:

$$\text{mg total dissolved solids/ 100 mL} = (A - B) \times 100 \times 1000 / \text{Volume of sample (mL)}$$

Where,

$$A - \text{Weight of the dried residue + dish (g)}$$

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Total Dissolved Solids (TDS) (Simplified Procedure):

This test should be carried out to estimate amount of total dissolved solids available in the effluent.

S. No.	Description	Details
01	Purpose of test	To estimate amount of total dissolved solids in the effluent.
02	Target value	< 350 mg /100 ml
03	Equipments required	Portable TDS meter
04	Other equipments	Beaker 100-150 ml
05	Quantity of sample	100 ml
06	Frequency of sampling	90 days
07	Testing spot	Railway laboratory
08	Staff required	01

Procedure:

1. Ensure availability of all the equipments and effluent in the lab
2. Fill the sample effluent in the beaker and mix it properly.
3. Keep the TDS meter ON and put electrode in the beaker up to marking line.
4. TDS value will appear on the display of TDS meter. It should be noted.
5. Repeat the test again and record the TDS value in the register.

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4. Total Volatile Solids (TVS):

This test should be carried out to estimate amount of total volatile solids available in the effluent.

In addition to equipments required in TDS test, muffle furnace is also used in this test

MUFFLE FURNACE

S. No.	Description	Details
1.	Purpose of test	To estimate amount of total volatile solids available in the effluent.
2.	Target value	< 500 mg /100 ml
3.	Equipments required	Electronic weighing balance, pipettes, Silica crucibles, Hot air oven, muffle furnace, desiccators, filter assembly.
4.	Consumable	Self indicating silica gel. Whatman Glass wool filters
5.	Quantity of sample	25 ml
6.	Electricity requirement	Yes
7.	Frequency of sampling	90 days
8.	Testing spot	Railway laboratory
9.	Staff required	1

Procedure:

1. Heat an empty and clean silica crucible at 550°C for 1/2 hour in muffle furnace.
2. Cool it in desiccator to room temperature and take the initial weight.
3. Pipette a measured volume of well mixed sample (25 ml) to a pre-weighed silica crucible.

4. Keep the silica crucible in a hot air oven at 103 – 105°C; keep till the water gets dried.
5. Remove the silica crucible and keep it in muffle furnace at 550°C for 1 hour.
6. Keep the silica crucibles in desiccators.
7. Note the final weight

Calculations: $\text{mg total volatile solids} / 100 \text{ mL} = (A - B) \times 100 / \text{Volume of sample (mL)}$

Where,

A – Total solids (mg)

B - Weight of the dried residue + dish (g) - Weight of dish (g) then convert the value to milli gram

Note : The total solids amount is calculated as described above.

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5. Chemical Oxygen Demand (COD) Test:

This test should be carried out to estimate COD of the effluent to ensure environmental parameters.

S. No.	Description	Details
1.	Purpose of test	To estimate COD of the effluent to ensure environmental parameters.
2.	Target value	< 2000 mg O ₂ / Litre
3.	Equipments	COD digester along with digestion tubes and condensers, Burettes, magnetic stirrer, magnetic bars.
4.	Consumable required	Sulphuric Acid, Silver sulphate, mercuric chloride, potassium dichromate, ferroin indicator, ferrous ammonium sulphate.
5.	Quantity of sample	5-10 ml of sample
6.	Electricity requirement	Yes
7.	Frequency of sampling	90 days
8.	Testing spot	Govt. approved Labs./DRDE
9.	Staff required	01

Reagents:

- A. Standard potassium dichromate solutions (0.04167M)** – Dissolve 12.259 g K₂Cr₂O₇, previously dried at 150°C for 2 hours, in distilled water and dilute to 1000 mL.
- B. Sulfuric acid reagent:** Add Ag₂SO₄ crystals or powder to concentrate H₂SO₄ at the rate of 5.5 g Ag₂SO₄/ Kg of H₂SO₄. Let it stand 1 or 2 days to dissolve.
- C. Ferroin indicator solution :** Dissolve 1.485 g of 1, 10 – phenanthroline monohydrate and 695 mg FeSO₄.7 H₂O in distilled water dilute to 1000 mL.
- D. Standard Ferrous ammonium Sulphate titrant (0.25)** – Dissolve 98 g of Fe(NH₄)₂.6 H₂O in distilled water. Add 20 mL con. H₂SO₄, cool and dilute to 1 litre standardize this solution on the day of estimation against Standard potassium dichromate solution.
- E. Mercuric sulphate - HgSO₄** crystals or powder.
- F. Sulfamic acid** – interference of nitrates is to be eliminated.

G. Potassium Hydrogen Phthalate: lightly crush and then dry KHP to constant weight at 110°C. Dissolve 425 mg in distilled water dilute to 1000 mL. KHP has a theoretical COD of 1.176 mg O₂ / mg and this solution has a theoretical COD of 500µg O₂ / mL.

Procedure:

- Take 5 ml of sample and dilute it to 50 ml. (Note: In case higher COD samples dilute 100 times).
- Add 100 mg HgSO₄, several glass beads/ Chemstones, and very slowly add 5.0 ml of sulfuric acid, with mixing to dissolve HgSO₄.
- Cool while mixing to avoid possible loss of volatile materials.
- Add 25.0 mL 0.04167 M K₂Cr₂O₇ solution and mix.
- Attach condenser to digestion tubes and cool the contents by swirling and mixing in running tap water.
- Add sulfuric acid reagent (20 mL) through open end of condenser.
- Continue swirling and mixing while adding sulfuric acid reagent.
- Cover open end of condenser with a small beaker / aluminum foil to prevent foreign material from entering refluxing mixture and to avoid possible loss of volatile materials.
- Do reflux / digestion for 2 hours @ 150°C.
- Cool and wash down condenser with distilled water.
- Disconnect reflux condenser and dilute mixture to about twice its volume with distilled water. Cool to room temperature.
- Titrate excess K₂Cr₂O₇ with FAS, using 0.10 – 0.15 mL (2 – 3 drops) ferroin indicator.
- Observe sharp colour change from blue green to reddish brown that persists for 1 min or longer.
- Duplicate determinations should agree within 5% of their average.
- In the same manner reflux and titrate a blank containing the reagents and a volume of distilled water equal to that of sample.

Determination of a STD solution:

Evaluate the technique and quality of the reagents by conducting the test on a standard KHP solution.

Calculations:

$\text{COD (mg/ml)} = (A-B) \times N \times 8000 \times \text{Dilution Factor} / \text{mL of sample}$

Where,

A – mL of FAS used for blank

B – mL of FAS used for sample

N – Molarity of FAS (0.25)

8000 – Mill equivalent weight of Oxygen.

Precautions:

1. Mix reflex mixture thoroughly before applying heat to prevent local heating of flask bottom and a possible blowout of flask contents.
2. Make dilution according to the COD range of the samples.
3. Ensure complete dissolution of silver sulphate in sulphuric acid reagent.
4. Store sulphuric acid reagent in a amber colored bottle in dark or wrap it with aluminum foil.
5. Wear lab coat and hand gloves during handling of acids.
6. Always add acid along the wall of the container to the water.
7. Do titration once the content reaches room temperature.

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6. Fecal Coli Forms count:

This test should be carried out to estimate the fecal coli form bacteria count on bio toilet effluent.

Laminar air flow chamber

Incubator

Micropipette

S. N	Description	Details
1.	Purpose of test	To estimate the fecal coli form bacteria count on effluent.
2.	Target value	> 99% reduction (Less than 10^8 /100 ML)
3.	Equipments	Laminar Air Flow Chamber, Incubator
4.	Consumables	Test tubes, FC media plates, spreaders, pipettes, conical flasks, Glass marker, magnetic stirrer.
5.	Quantity of sample	10 ml
6.	Electricity requirement	Yes
7.	Frequency of sampling	90 days
8.	Testing spot	Only in Govt. Labs.
9.	Staff required	01

Procedure :

Media & Plate Preparation:

A. Media Composition :

S. No.	Constituents	Amount (g/ L)
1.	Lactose	12.5
2.	Tryptone	1
3.	Protease Peptone	5
4.	Sodium chloride	5
5.	Yeast extract	3
6.	Bile salts	1.5
7.	Aniline blue	0.1
8.	Distilled water	1000 ml

1. Boil the media for 30 minutes, do not sterilize the media.
2. Add approximately 25 ml of cooled FC medium to each Petri dish. Upon solidification of the media close the lids and keep it for further use.

B. Dilution preparation:

1. Label 90 ml water blank conical flask as **No.1**(10^{-1}) and subsequent 9 ml water blank test tubes as **No. 2** (10^{-2}), **No. 3** (10^{-3}), **No. 4** (10^{-4}), **No. 5** (10^{-5}), **No. 6** (10^{-6}), **No. 7** (10^{-7}), **No. 8** (10^{-8}), with a glass marker.
2. Take 10 ml of well mixed effluent sample and add to 90 ml water blank to make (10^{-1}) dilution. (Note: Use sterile pipettes for transfer and carry out dilution in laminar air flow chamber)
3. Vigorously shake the dilution in a magnetic stirrer for 2 minutes to obtain uniform suspension.
4. Transfer 1 ml of suspension from No.1 into No.2 tube to make (10^{-2}) dilution and shake it vigorously for 1 min.
5. Make further dilutions as prepared above.

Note: Do estimation of fecal coli form count in human excreta for comparison and calculation.

C. Inoculation:

1. Place the spreader into a beaker containing 95% ethyl alcohol. Remove the spreader and pass it through burner flame. Burn off the alcohol completely and cool the rod for 30 to 45 seconds.
2. Transfer 0.1 ml (100 μ l) of suspension each from different dilution into 3 plates. Place it in the center of the plate.
3. Remove the Petri dish cover, with one hand touch the spreader gently on the surface of the agar and move it back and forth while rotating the plate with the other hand.
4. When the suspension gets dried or absorbed on the agar medium, replace the cover.
5. Incubate the plates in an upright position in an incubator at 44.5°C for 48 hours.

Observations :

1. Observe the plates for number and distribution of Blue/ Bluish tinged colonies. Do not consider white or creamy or dull white colonies.
2. Select the plate from the appropriate dilution which contains colonies in the range of 30 to 300 and make the plate counts.
3. Determine the average of the triplicate colony count.
4. Record your results.

Calculation:

Calculate fecal form count in terms of number of CFU / 100 ml by applying the following formulas.

$$\text{Fecal coli form count (CFU / 100 ml)} = \frac{\text{Mean plate count} \times \text{Dilution factor} \times 100}{\text{Volume of sample}}$$

$$\% \text{ reduction} = \frac{(\text{count in human excreta} - \text{count in sample}) \times 100}{\text{count in human excreta}}$$

Example

If 160 colonies were counted (average of 3 replicate) in 10^{-3} dilution the fecal coli form count would be:

$$\begin{aligned} \text{FC count (CFU / 100 ml)} &= \frac{160 \times 1000 \times 100}{0.1} \\ &= 1.6 \times 10^8 \end{aligned}$$

Precautions:

1. The sample should be mixed homogenously.
2. Each dilution must be thoroughly shaken before removing an aliquot for subsequent dilution.
3. Use separate sterile pipettes or tips for each dilution.
4. All dilution making and inoculation should be done in laminar air flow chamber.
5. The plates should be incubated in a upright position.
6. 25-30 ml media should be added to avoid cracking of medium during incubation.

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(B) Testing Scheme for Bacteria Culture (Microbial Inoculum).

Following parameters have been recommended for testing of Microbial Inoculum produced in Inoculum generation Plant .

1. pH Value Test:

This test should be carried out to measure pH value of the inoculum.

S. No.	Description	Details
1.	Purpose of test	To measure pH value of the inoculum.
2.	Target value	6.5 – 7.5. pH
3.	Equipments required	Table top pH meter / pH indicator strips, magnetic stirrer.
4.	Consumables	pH calibration buffer (4.0, 7.0, 10.0), magnetic stirrer bars.
5.	Quantity of sample	50 – 100 ml.
6.	Electricity requirement	yes
7.	Frequency of sampling	Once in a week.
8.	Testing spot	Railway laboratory
9.	Staff required	01

Procedure:

- Take 50 – 100 ml of well mixed inoculum in a beaker.
- Put a magnetic bar and keep the beaker on a magnetic stirrer and mix it continuously.
- Wash the electrode and Temperature compensation rod with distilled water and wipe it with tissue paper.
- Put the electrode and automatic temperature compensation rod into the sample and keep it until stable reading appears in the display; note the reading.
- Discard the sample and wash the electrode and automatic temperature compensation rod with distilled water and wipe it with tissue paper.
- Keep the electrode back in the container and put them in stand.

Precaution: Do calibration with appropriate buffers before taking readings.

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2. Bio-gas Test:

This test should be carried out to measure amount of Biogas produced during fermentation in the Inoculum generation plant.

S. No.	Description	Details
1.	Purpose of test	To measure amount of biogas produced during fermentation.
2.	Target value	-
3.	Equipments	Gas flow meter / water displacement assembly.
4.	Consumables	Glass wares, cork, Glass 'L' bends
5.	Quantity of sample	Not applicable
6.	Electricity requirement	No.
7.	Frequency of sampling	Daily
8.	Testing spot	At Plant/Site
9.	Staff required	01

Procedure:

- Take 2000 ml of well mixed inoculum in a 3 L beaker.
- Set the water displacement assembly.
- Incubate the set-up in BOD incubator @ 30 °C for 24 hours.
- Measure the volume of water displaced.
- The amount of water displaced is equal to amount of gas produced.

Precaution: Make sure water column in the exit tube of the water beaker.

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3. Percentage Methane Test:

This test should be carried out to measure the methane content of biogas produced during fermentation in the Inoculum generation plant.

S. N	Description	Details
1.	Purpose of test	To measure the methane content of biogas produced during fermentation.
2.	Target value	40 – 70 %
3.	Equipments required	Gas chromatograph / online gas analyzer.
4.	Consumables	Serum bottles, Aluminum crimp, Butyl rubber septum, Gas tight syringe (50 – 250 µl), Nitrogen cylinder, Zero air/O ₂ cylinder, Hydrogen cylinder, standard methane gas.
5.	Quantity of sample	Not applicable
6.	Electricity requirement	Yes.
7.	Frequency of sampling	weekly
8.	Testing spot	Railway laboratory
9.	Staff required	01

Procedure:

- Take 2000 ml of well mixed inoculum in a 3 L beaker.
- Set the water displacement assembly.
- Incubate the set-up in BOD incubator @ 30 °C for 24 hours.
- Fill a 30 ml serum vial with water, put a butyl rubber septum, aluminum crimp and crimp it all the side.
- Inject the biogas formed in the water displacement assembly into the serum vial through a silicone tube fitted with needle.
- Allow the exit of water through another needle inserted in the butyl rubber septum.
- Collect at least 15 ml of biogas.
- Store it in a inverted position until analysis.
- *Note: Make sure, gas bubble should not be there in the serum bottle before biogas collection.*
- Inject 100 µl standard methane gas and record the elution time of the gas.
(*Note: Elution time is the time elapsed between injection of sample and response of the peak at top*).
- Inject 100 µl of the biogas sample. Record the % methane in biogas.

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4. MPN Count for Methanogens:

This test should be carried out to estimate the methane producing bacteria count in inoculum.

S. No.	Description	Details
1.	Purpose of test	To estimate the methane producing bacteria count in inoculum.
2.	Target value	>1000 / ml
3.	Equipments required	Gas chromatograph, Gassing manifold, hotplate stirrer, laminar air flow chamber, incubator, micropipettes, hand operated crimper, hand operated Decapper.
4.	Consumables	Serum bottles, Aluminum crimp, Butyl rubber septum, Gas tight syringe (50 – 250 µl), Nitrogen cylinder, CO ₂ cylinder, Hydrogen cylinder, sterile Disposable syringes, Resazurin and chemicals listed below in the table.
5.	Quantity of sample	~ 50 ml
6.	Electricity requirement	Yes.
7.	Frequency of sampling	Monthly
8.	Testing spot	In Govt. approved Laboratory Only
9.	Staff required	02

Procedure:

Growth medium preparation:

Sl No	Constituents	Amount (g / L)
1	KH ₂ PO ₄	0.3
2	K ₂ HPO ₄	0.3
3	(NH ₄) ₂ SO ₄	0.3
4	NaCl	0.6
5	MgSO ₄ .7H ₂ O	0.13
6	CaCl ₂ .2H ₂ O	0.008
7	FeSO ₄	0.002
8	Yeast extract	2.0
9	Trypticase	2.0
10	Trace element solution	1.0 ml
11	Trace vitamin solution	1.0 ml
12	Sodium acetate	2.5
13	Sodium formate	2.5
14	Resazurin	1.0 ml
15	Reducing agent	12.5 ml
16	pH	6.8
17	Gas phase (N ₂ :H ₂)	80:20

Trace Element Solution:

SN	Constituents	Amount (g / 100 ml)
1	Nitrilo triacetic acid	1.5
2	MgSO ₄ .7H ₂ O	3.0
3	MnSO ₄ .2H ₂ O	0.5
4	NaCl	1.0
5	CoCl ₂	0.1
6	FeSO ₄ .7H ₂ O	0.1
7	CaCl ₂ .2H ₂ O	0.1
8	ZnSO ₄	0.1
9	CuSO ₄ .5H ₂ O	0.01
10	AlK(SO ₄) ₂	0.01
11	H ₃ BO ₃	0.01
12	Na ₂ MoO ₄ .2H ₂ O	0.1

Note: Dissolve nitrilotriacetic acid with KOH to pH 6.5 and then add other salts.

Trace Vitamin Solution:

Sl No	Constituents	Amount (mg / 100 ml)
1	Biotin	2.0
2	Folic acid	2.0
3	Vitamin B ₁₂	0.1
4	Pyridoxine HCl	10.0
5	Thiamine	5.0
6	Riboflavin	5.0
7	Nicotinic acid	5.0
8	D,L- Calcium pantothenate	5.0
9	p- aminobenzoic acid	5.0
10	Lipoic acid	5.0

Reducing Agent:

Sl No	Constituents	Amount (g / 100 ml)
1	Cysteine hydrochloride	2.0
2	Sodium sulphide	2.0

Adjust pH of Cysteine hydrochloride to 9.0

Indicator Solution:

Sl No	Constituents	Amount (g / 100 ml)
1	Resazurin	0.10

1. All the ingredients of the medium were added in the required amount except the heat labile components and half the volume (6.25 ml) of reducing agent.
2. Heat the medium on a hot plate or heating mantle at 70 °C and sparge the medium continuously with O₂ free N₂ till the medium becomes colorless.
3. Take 60 ml serum vials and replace the air with nitrogen for 2 min.

4. Dispense media (22.5 ml) in individual vials and seal it with butyl rubber and aluminum crimps and seal it tightly.
5. Do autoclaving for 15 min at 121 °C for 20 min.
6. Let the medium to cool to room temperature and add rest of the amount of reducing agent (6.25 ml) with syringe and heat the butyl rubber seal for 10 sec with a spirit lamp.

Note: *The medium should be colourless.*

Dilution Preparation:

1. Label 22.5 ml anaerobic diluent medium containing serum vials as **No. 1** (10^{-1}) **No. 2** (10^{-2}), **No. 3** (10^{-3}), **No. 4** (10^{-4}), **No. 5** (10^{-5}), **No. 6** (10^{-6}), **No. 7** (10^{-7}), **No. 8** (10^{-8}) with a glass marker.
2. Take 2.5 ml of well mixed effluent sample and add to 22.5 ml water blank to make 10^{-1} dilution. Vigorously shake the dilution for 15 seconds to obtain uniform suspension.
3. Transfer 2.5 ml of suspension from No. 1 into No. 2 tube to make 10^{-2} dilution with a sterile syringe and shake it vigorously for 15 seconds.
4. Make further dilutions as prepared above.

Inoculation:

- For each dilution prepare 5 replicate serum bottles for inoculation. Total 45 medium bottles are needed.
- Transfer 2.5 ml of suspension each from different dilutions into set of 5 serum vials with medium by sterile syringe and heat the butyl rubber septum immediately.
- Incubate the serum vials in an upright position along with uninoculated control in an incubator at 30 °C for 20 days.

Observations:

- Analyze the gas from the head phase for presence of methane as described above with gas chromatograph.
- Note the number of positive and negative vials for each dilution.
- Methanogenic MPN is computed on the basis of bottles showing positive test. Dilution factor is taken into consideration while calculating the MPN.
- Record your results

Calculation:

1. To calculate MPN of methanogens in the original sample
 - Select as P1: the number of positive tubes in the least concentrated dilution in which the greatest number of tubes is positive; and p2 and p3 represent the numbers of positive tubes in the next two higher dilutions.
 - Use table of MPN for use with 5 tubes per dilution, then find the row of numbers in MPN table in which p1 and p2 correspond to the values observed experimentally. Follow that row of numbers across the table to the column headed by the observed value of p3. The figure at the point of intersection is the most probable number of organisms in the quantity of the original sample represented in the inoculum added

in the second dilution (p2). Multiply this figure by the appropriate dilution factor to obtain the MPN for the original sample.

Example

Suppose the following observations were made:

Dilution	No. of positive
Neat	5
10^{-1}	5
10^{-2}	5
10^{-3}	5
10^{-4}	3
10^{-5}	1
10^{-6}	0

In this series, $p_1 = 5$, $p_2 = 3$, $p_3 = 1$. For this combination of p_1 , p_2 and p_3 , the MPN table gives 1.1 as the most probable number of organisms in the quantity of the inoculum applied in the 10^{-4} (p_2) dilution. Multiplying this number with dilution factor 10^4 gives 1.1×10^4 as the MPN for the original sample (2.5 ml). Convert the value to 1 ml by dividing with 2.5 i.e. $1.1 \times 10^4 / 2.5 = 4.4 \times 10^3$

Precautions:

1. The sample should be mixed homogeneously
2. Each dilution must be thoroughly shaken before removing an aliquot for subsequent dilution
3. Use separate sterile pipettes or tips for each dilution
4. Use sterile pipettes for transfer,
5. Carry out dilution as quick as possible and immediately seal the butyl rubber septum and aluminum crimps. Alternatively the dilution can be done in anaerobic glove box.
6. The Heat labile compounds should be added after autoclaving.
7. Trace element solution and reducing agent solution should be prepared separately and should be autoclaved.
8. Trace vitamin solution should be sterilized by filtration.
9. At least three consecutive dilutions should be inoculated

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Category 0: Unacceptable tube combinations, having, under normal conditions, the least chance to be obtained. Combinations not mentioned in the table also belong to this category.

Category 1: Acceptable combinations, having, under normal conditions, the greatest chance to be obtained. Category 1 combinations should normally be obtained in 95% of the cases.

Category 2: Combinations that may be accepted if decisions based on them are not of utmost importance. They are less likely to be obtained than category 1 combinations, but not as unlikely as category 0. In 99% of the cases a combination belonging to category 1 or 2 should be obtained

Date of Sampling:

Date of Lab Testing:

Coach No.	Variant of Bio-digester	Observations/Results					
		pH Value	Total Solids (mg/100ml)	TDS (mg/ 100ml)	COD Level (ppm)	Volatile Solids (mg/100ml)	Faecal Coliform Count (MPN/ 100ml)

Signature of Evaluating Authority :
 Name & Designation :
 Seal :

Bio-Toilet Maintenance Proforma

Lavatories arrived choked and cleared during Maintenance						
S.No.	Date	Coach No. Manufacturer/ CLASS/ VARIANT	Lavatory Number	Reason for chocking	Signature of Railway Engineer	Signature of representative of AMOC Agency
1.						
2.						
3.						
4.						
5.						
6.						
7.						
8.						
9.						
10.						
11.						
12.						
13.						
14.						
15.						
16.						
17.						
18.						
19.						
20.						
21.						
22.						

Remarks: Description of choked lavatory & leakages.

Field Trial of Bio-Toilet (IR-DRDO Technology) on Indian Railway Passenger Coaches

Train No..... Rake No..... Coach No.Date.....Name of SE/ JE R/ Maint.....

S.N	Parameters to be checked	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
1.	Thorough inspection of mounting / securing arrangement for the bio-digester tank.																
2.	Any leakage in joints/ connections in the complete system including retention tank, water pipe line.																
3.	Adequacy of the provisions made for segregation of non bio-degradable items.																
4.	Functionality of flush buttons/ lever & water tabs																
5.	Adequacy of the flushing of the pan																
	(a) By using pressurized flushing (b) By normal manual flushing																
6.	Any mal odour/ stench: (yes/No), if yes, please specify Light, medium or Heavy																
7.	Overall cleanliness level of the toilet room.																
8.	Notices for users and maintenance personal in Hindi/English & regional language of the originating & destination station.																
9.	Emergency operation of flush arrangement without power / air supply.																
10.	Performance of non-bio-degradable waste ejection system its effectiveness and reliability.																
	a. Working satisfactorily (Yes/ No), if no, give reasons and remarks for malfunctioning. b. Choking of toilet pan (Yes/ No), if yes give reasons and remarks for choking.																
11.	Tank evacuation date (if any) with coach No. and date.																
12.	Details of attentions / maintenance required with coach no. , date and time taken for maintenance.																
13.	Adequacy of strength & capacity of waste treatment tank																
14.	Choking of P-Trap (Yes/ No), if yes, give reasons and it has been attended or not.																
15.	Check operation of flapper / Ball valve																

Provision of Stickers in Bio-Toilets

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List of Govt. Approved Environmental Laboratories

ENVIRONMENTAL LABORATORIES WITH VALID RECOGNITION UNDER SECTION 12(1) B OF THE ENVIRONMENT (PROTECTION) ACT, 1986

The Environmental laboratories (Govt./Semi-Govt./Public Sector Undertakings/Educational Institutes) with valid recognition (updated upto 15 April, 2012) under the Environment (Protection) Act, 1986 are as below :

S. No.	Name of laboratory	Group of parameters	Gazette notification no. and date	Validity upto
1.	Goa State Pollution Control Board EDC Plaza Patto, Panaji, Goa-403001	Physical, Chemical, Microbiological, and Air Pollution parameters for analysis of ambient air, source emission and micro-meteorological parameters.	Legal 42(3)/87 dated 11 th October, 2007	30 th September, 2012
2.	Centre Laboratory Madhya Pradesh Pollution Control Board "Paryavaran Parisar", E-5, Arera Colony Bhopal-462 018. M. P.	Physical. Chemical, Microbiological. Toxicological and Air Pollution parameters for analysis of ambient air, sources emission and micro-meteorological parameters.	Legal 42(3)/87 dated 1 st October, 2007	30 th September, 2012
3.	Eco-Auditing Laboratory National Botanical Research Institute Pana Pratap Marg, Lucknow-228 001 U.P.	Physical. Chemical, Microbiological, and Air Pollution parameters for analysis of ambient air, source emission and micro-meteorological parameters.	Legal 42(3)/87 dated 1 st October, 2007	30 th September, 2012
4.	Zonal Office Laboratory, Central Pollution Control Board Synergy House-II Gorwa Subhanpura Road Subhanpura Vadodara-390 023 ,Gujarat .	Physical, Chemical, Microbiological. Toxicological and Air Pollution parameters for analysis of ambient air, source emission and micro-meteorological parameters.	Legal 42(3)/87 dated 1 st October, 2007	30 th September, 2012
5.	Central Laboratory Pollution Control Board. Assam Bamunimaidan Guwahati-781 021. Assam	Physical. Chemical. Microbiological, Toxicological and Air Pollution parameters for analysis of ambient air, source emission and micro-meteorological parameter.	Legal 42(3)/87 dated 1-February, 2008	31 st January, 2013
6.	Environmental Management Division Central Pulp & Paper Research Institute, Post Box No. 174, Paper Mills Road, Himmat Nagar Saharanpur 247001, U.P.	Physical, Chemical. Microbiological, Toxicological. biological and Air Pollution parameters for analysis of ambient air. source emission and micro-meteorological parameters.	Legal 42(3)/87 dated 10 th February, 2008	31 st January, 2013

S. No.	Name of laboratory	Group of parameters	Gazette notification no. and date	Validity upto
7.	Central Laboratory National Fertilizers Limited Vijaipur Unit, Vijaipur-743111 Tehsil-Raghogarh Dist.Guna, M. P.	Physical. Chemical, Microbiological. Toxicological and Air Pollution parameters for analysis of ambient air. source emission and micro- meteorological parameters.	Legal 42(3)187 dated	31 st January, 2013
			1 st February, 2008	
8.	Central Laboratory Andhra Pradesh Pollution Control Board A-3. Industrial Estate Sanathnagar Hyderabad-500 018 Andhra Pradesh	Physical. Chemical. Organics Microbiological. Toxicological. Biological. Hazardous Waste and 1Jr Pollution parameters for analysis of ambient air. Source emission and micro-meteorological parameters.	Legal 42(3)/87 dated 1 st February, 2008	31 st January, 2013
9.	Gujarat Pollution Control Board, Regional Office -Surat 338, Typical First Floor Belgium Square Sliver Plaza Complex, Opp. Linear Bus Stand, Ring Road Surat-395 003, Gujarat	Physical. Chemical. Organics Microbiological. Toxicological. and Air Pollution parameters for analysis of ambient air. source emission. noise. vehicular emission and micro-meteorological parameters	Legal 42(3)187 dated 1 st February. 2008	31 st January. 2013
10.	Gujarat Pollution Control Board Regional Office Laboratory- Vadodara Geri Compound Race Course. Opp. S. T. Depot Vadodara-390 007 Gujarat	Physical. Chemical, Organics Micro biological. Toxicological. and Air Pollution parameters for analysis of ambient air, .source emission. noise. vehicular emission and micro-meteorological parameters	Legal 42(3)187 dated 1 st February, 2008	31 st January, 2013
11.	Regional Office Laboratory- Rajkot Gujarat Pollution Control Board Race Course. Ring Road Near Union Bank, Rajkot-360 001 Gujarat	Physical. Chemical. Organics Microbiological. Toxicological. and Air Pollution parameters for analysis of ambient air. source emission. noise. vehicular emission and micro-meteorological parameters	Legal 42(3)/87 dated 1 st February. 2008	31 st January. 2013
12.	Environmental Laboratory Central Mine Planning & Design Institute limited Gondwana Place ,Kanke Road, Ranchi-834 008. Jharkhand	Physical. Chemical, Microbiological, and Air Pollution parameters for analysis of ambient air. source emission and micro-meteorological parameters.	Legal 42(3)/87 dated 1 st February. 2008	31 st January. 2013
13.	Pollution Control Research institute, Bharat Heavy Electricals Limited Ranipur. Haridwar-249 403 Uttarakhand	Physical. Chemical. Microbiological. Toxicological and Air Pollution parameters for analysis of ambient air. source emission. and, micro- meteorological parameters.	Legal 42(3)/87 dated 1 st April. 2008	31 st March. 2013

S.No	Name of laboratory	Group of parameters	Gazette notification no. and date	Validity upto
14.	Central Laboratory Central Pollution Control Board Parivesh Bhawan East Arjun Nagar Delhi-110032	Physical. General and non- Metallic. Metals. Organics. Microbiological. Toxicological Biological, Hazardous Wastes, Soil, Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission. noise and micro-meteorological parameters,	Legal 42(3)/87 dated 1 st October, 2009	30 th September. 2014
15.	Central Laboratory Punjab Pollution Control Board Vatavaran Bhawan. Patiala, Punjab-147001	Physical, General Chemical and non-Metallic, Metals, Organics, Microbiological, Toxicological Biological, Hazardous wastes. Soil. Sludge, Sediments and Air Pollution Parameters for analysis of ambient air, source emission, noise and micro-meteorological & vehicular emission monitoring parameters,	Legal 42(3)/87 dated 15 th January. 2010	14 th January. 2015
16.	Regional Laboratory Maharashtra Pollution Board Control Jog Centre, 3rd Floor, Puna-Mumbai Road, Shivaji Nagar. Pune-411003	Physical, General, Chemical and non-Metallic, Metals, Organics, Microbiological, Toxicological, Hazardous wastes, Soil, Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological parameters.	Legal 42(3)/87 dated 15 January, 2010	14 th January , 2015
17.	Zonal Laboratory Central Pollution Control Board, Zonal office , Kolkata 502, South end Conclave, 1582 Rajdanga Main Road Kolkata 700107	Physical, General, Chemical and non-Metallic. Metals, Organics, Microbiological, Toxicological Biological, Hazardous Wastes, Soil, Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological & vehicular emission monitoring parameters,	Legal 42(3)/87 dated 15 January, 2010	14 th January. 2015
18.	Environment Protection Training and Research Institute (EPTRI), 91/4, Gachi Bowli. Hyderabad. 500032 Andhra Pradesh	Physical, General, Chemical and non-Metallic, Metals, Organics, Microbiological. Toxicological Biological, Hazardous Wastes, Soil, Sludge. Sediments and-Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological Parameters,	Ler-142(3)/87 dated 20 September. 2010	19 th September, 2015
19.	P.G. Department of Environment Management Chhatrapati Shahu Institute of Business Education and Research (SIBER). University Road, Kolhapur. 416004 Maharashtra	Physical, General, Chemical and non-Metallic, Metal, Organics, Microbiological, Toxicological, Biological, and Air Pollution parameters for analysis of ambient air, source emission, and micro-meteorological parameters.	Legal 42(3)/87 dated 20 September. 2010	19 th September. 2015

S. No.	Name of laboratory	Group of parameter's	Gazette notification no. and date	Validity upto
20.	Punjab Bio-Technology Incubator Agri. And Food. testing Laboratory SCO: 7 & 8 (Top Floor). Phase.V, SAS Nagar. Mohali-160059 Punjab.	Physical, General. Chemical and non-Metallic. Metals. Organics, Microbiological. Toxicological. Hazardous Wastes, Soil. Sludge. Sediments and Air Pollution parameters for analysis of ambient air, source emission, and micro-meteorological parameters.	Legal 42(3)187 dated 20 th September. 2010	19 th September, 2015
21.	Regional Laboratory Maharashtra State Pollution Control Board 6th Floor, "Udyog Bhawan" Civil Lines Nagpur-440001 Maharashtra	Physical. General. Chemical and non-Metallic, Metals. Organics, Microbiological, Toxicological Biological. Hazardous Wastes, Soil, Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission. noise and micro-meteorological parameters.	Legal 42(3)/87 dated 2 nd January. 2011	26 th January, 2016
22.	Regional Laboratory Maharashtra State Pollution Control Board 1 st Floor Udyog Bhawan. Rathi Chowk, Trimbak Road, Nashik-422007 Maharashtra	Physical, General, Chemical and non-Metallic, Metals, Organics, Microbiological, Toxicological Biological, Hazardous wastes, Soil, Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological parameters.	Legal 42(3)187 dated 27 th January, 2011	26 th January, 2016
23.	Regional Laboratory Maharashtra State Pollution Control Board "Paryavaran Bhavan" A-4/1, Chikalthana MIDC, Behind Dhoot Hospital, Aurangabad-431210 Maharashtra	Physical, General, Chemical and non-Metallic, Metals, Organics, Microbiological. Toxicological Biological, Hazardous Wastes, Soil. Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological Parameters.	Legal 42(3)187 dated 27 th January. 2011	26 th January. 2016
24.	Central Laboratory Uttar Pradesh Pollution Control Board PICUP Bhawan, 3 rd Floor, B-Block Vibhuti Khand, Gomati Nagar Lucknow.226010 U.P.	Physical, General, Chemical and non-Metallic, Metals, Organics. Microbiological. Toxicological. Soil, Sludge. Sediments and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological parameters.	Legal 42(3)187 dated 31 st January, 2012	•
25.	Regional Laboratory M.P. Pollution Control Board Plot No, 455/456, Vijay Nagar. Jabalpur-482002 Madhya Pradesh	Physical. General, Chemical and non-Metallic, Metals, Organics, Microbiological, Toxicological, Soli. Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological Parameters.	Legal 42(3)/87 dated 31 st January. 2012	•

S. No.	Name of laboratory	Group of parameters	Gazette notification no. and date	Validity upto
26.	Regional Laboratory M.P. Pollution Control Board Sc-17, Bharatpuri Ujjain-456010 Madhya Pradesh	Physical, General, Chemical and non-Metallic, Metals, Organics, Microbiological, Toxicological, Soil, Sludge, Sediments and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological parameters.	Legal 42(3)187 dated 31 st January, 2012	•
27.	Regional Laboratory M.P. Pollution Control Board Scheme No.-78 Part-II. Aranya Nagar Indore 452010 Madhya Pradesh	Physical, General, Chemical and non-Metallic, Metals, Organics, Microbiological. Toxicological. Soil. Sludge, Sediment; -and Air Pollution parameters for analysis of ambient air, source emission, noise and micro-meteorological parameters.	Legal 42(3)/87 dated 31 st January, 2012	•

- Environmental laboratories and the Government Analysis so mentioned shall remain valid only up to 31st Aug 2012 as per provision of office memorandum dated 12th AUG 2011 Issued by ministry of Environment & forest to acquire the required accreditation / certification to achieve eligibility under the essential pre-requisite or office memorandum.