

# A HOMOGENEOUS FLUORESCENCE POLARIZATION ASSAY FOR NEUROKININ-1 RECEPTORS USING Cy3B LABELLED SUBSTANCE P

Alison J. Harris, Sarah L. Cox, Malcolm Allen, Julian Long, \*Christopher G. Norey, and D. Dougal Burns.

Amersham Biosciences UK Limited, Amersham Place, Little Chalfont, Buckinghamshire HP7 9NA England. Telephone: +44 (0)29 2052 6439, Fax: +44 (0)29 2052 6230, e-mail: christopher.norey@eu.amershambiosciences.com

## Introduction

Fluorescence Polarization (FP) is a homogeneous technique that is now widely used for high throughput screening (HTS). We have labelled Substance P (Sub P) with Cy<sup>TM</sup>3B at the N-terminus or at the internal lysine position (see below). These peptide ligands together with FARCyte<sup>TM</sup> fluorescent plate reader have been used to develop 384-well and miniaturized 1536-well FP assays for the Neurokinin-1 (NK<sub>1</sub>) receptor.

Cy3B-(terminal)-Sub P: Cy3B-RPKPQQFFGLM  
Cy3B-(lysine)-Sub P: RPK(Cy3B)-PQFFGLM

## Method

Peptides were prepared by automated synthesis where acetylation at the N-terminus or ε-lysine facilitated site specific labelling. Cy3B was introduced via *in situ* activation and the labelled peptides purified by conventional reverse phase HPLC. The lyophilized products were characterized by UV spectroscopy and Maldi ToF mass spectroscopy. The affinities of Cy3B labelled Sub P ligands for the NK<sub>1</sub> receptor (Table 1) were determined by radioactive Scintillation Proximity Assay (SPA) and compared with the native ligand.

Table 1: Substance P ligand IC<sub>50</sub> values.

Compound	IC <sub>50</sub>
Sub P	0.17nM
Cy3B-(lysine)-Sub P	5nM
Cy3B-(terminal)-Sub P	2nM


CHO NK<sub>1</sub> membranes (5μg) were incubated in the presence of 70pM [<sup>35</sup>S]Sub P with varying concentrations of competitors; Sub P, Cy3B-(lysine)-Sub P or Cy3B-(terminal)-Sub P together with 0.5mg Type B PEI WGA SPA beads. Non-specific binding was determined in the presence of 25μM Sub P. Incubations were performed at room temperature for 3 hours in 50mM Tris, pH7.5 containing 0.5%(w/v) BSA, 1mM EDTA, 1mM MnCl<sub>2</sub>, 400μg/ml bacitracin, 20μg/ml leupeptin and 1mM Pefabloc<sup>TM</sup> SC. The assay plate was counted using the EG&G Wallac 1450 MicroBeta<sup>TM</sup> scintillation counter.

For 384-well FP assays, CHO NK<sub>1</sub> cell membranes (~1pmol/mg) were incubated with either of the Cy3B ligands (0.5nM) in the presence of varying amounts of Sub P, [Sar<sup>9</sup>Met(O<sub>2</sub>)<sup>11</sup>] Sub P or Neurokinin-A (NKA) for 3.5 hours at room temperature (20 – 25°C) in 50mM Tris buffer pH7.5 containing 0.5%(w/v) BSA, 1mM EDTA, 1mM MnCl<sub>2</sub>, 400μg/ml bacitracin, 20μg/ml leupeptin and 1mM Pefabloc SC. Incubations were performed in a total volume of 50μl in Corning black 384-well non-binding surface plates. Non-specific binding (NSB) was determined in the presence of 25μM Sub P. FP was then determined on FARCyte using a Cy3B optical filter/dichroic configuration. Specific polarization values were calculated by subtracting the polarization signal obtained in NSB wells (in the presence of 25μM Sub P) from the test well values. Further assays were performed as described with Cy3B-(terminal)-Sub P (0.5nM) only, in a total volume of 8μl in Greiner black 1536-well microtitre plates.

## Results

Maximal FP change in the competitive displacement assay signal was observed using Cy3B-(terminal)-Sub P but consistent results were obtained with both ligands (Figure 1). Similar performance using Cy3B-(terminal)-Sub P was also observed during further miniaturization of the assay to a 1536-well format (Figure 2).

Figure 1: NK<sub>1</sub> competitive binding assays (Sub P, [Sar<sup>9</sup>Met(O<sub>2</sub>)<sup>11</sup>] Sub P and NKA) in 384-well format.




CHO NK<sub>1</sub> cell membranes (30μg) were incubated with a) 0.5nM Cy3B-(terminal)-Sub P or b) 0.5nM Cy3B-(lysine)-Sub P in the presence of Sub P (■), [Sar<sup>9</sup>Met(O<sub>2</sub>)<sup>11</sup>]-Sub P (▲) or NKA (▼) as described in the text and FP was measured on FARCyte. Values are plotted as means of quadruplicates ± SEM.

Z' factor analysis, recently described by Zhang *et al*<sup>1</sup>, describes assays with a Z' factor between 0.5 and 1.0 to be reliable and robust for HTS. Here values of 0.68 and 0.51 were obtained for the 384-well and 1536-well formats respectively (Figure 3), indicating both are suitable for HTS purposes.


Table 2 compares both the 384-well and 1536-well assays for the NK<sub>1</sub> receptor. Miniaturization represented an approximate 6-fold reduction in membrane/ligand usage and assay volume. In addition, the ability to monitor competitive displacement and observe acceptable assay reliability was retained.

Figure 2: NK<sub>1</sub> competitive binding assay with Sub P in 1536-well format.



CHO NK<sub>1</sub> cell membranes (5μg) were incubated with 0.5nM Cy3B-(terminal)-Sub P and varying concentrations of Sub P as described in the text and FP was subsequently measured on FARCyte. Values are plotted as means of quadruplicates ± SEM.

Figure 3: Z' analysis for the NK<sub>1</sub> FP assays.



CHO NK<sub>1</sub> cell membranes (a) 30μg in 384-well plates and (b) 5μg in 1536-well plates were incubated with 0.5nM Cy3B-(terminal)-Sub P in the absence (■) or the presence (▲) of 25μM Sub P as described in the text and FP measured on 168 replicates of each on FARCyte.


In order to assess the applicability of such assays to HTS environments, the effect of DMSO was investigated. The assay tolerated up to 2.5% (v/v)

DMSO. At higher concentrations a 70% reduction in the binding assay signal was observed (Figure 4a). We conclude this to be an effect on binding efficacy in this assay as when DMSO was added after the completion of binding, signal was reduced by only 20% even in the presence of 10% DMSO (Figure 4b).

Table 2: Summary of FP assay miniaturization.

Assay conditions	384-well	1536-well	Reduction (fold)
Membrane protein	30μg	5μg	6
Cy3B-(terminal)-Sub P	25fmol	4fmol	6.25
Volume of the assay	50μl	8μl	6.25
Final Conc. of Cy3B-(terminal)-Sub P	0.5nM	0.5nM	-
Apparent IC <sub>50</sub> (Sub P)	4.3nM	6.5nM	-
Assay Z' factor	0.68	0.51	-

Figure 4: Effect of DMSO on the NK<sub>1</sub> FP assay.



CHO NK<sub>1</sub> membranes (30μg) were incubated with 0.5nM Cy3B-(terminal)-Sub P as described in the text together with varying concentrations (v/v) of DMSO (0.5, 1, 2.5, 5 or 10%) either (a) from the outset (■) or (b) following completion (▲) of the binding incubation. FP was subsequently measured on FARCyte. Values are plotted as means of quadruplicates ± SEM.

## CONCLUSIONS

- We have previously miniaturized a small molecule FP receptor ligand assay<sup>2</sup> and here have successfully miniaturized an FP receptor peptide ligand assay in 384 and 1536-well formats using a Cy3B label and FARCyte.
- The miniaturization significantly reduced reagent consumption and assay volume without compromising assay reliability.
- The apparent IC<sub>50</sub> values obtained in both formats were consistent and Z' values of >0.5 were also observed.
- These assays have used common receptor expression levels (~1pmol/mg) and low ligand concentrations (0.5nM), extending the potential of this technique to other small molecule or peptide ligand binding assays in HTS.

## References:

- Zhang, J.-H., *et al* (1999), *J. Biomol. Screening*, 2, 67 - 73.
- Harris, A.J., *et al* (2001), poster presentation at "Drug Discovery Technologies", Boston, MA,

Scintillation Proximity Assay (SPA) technology is covered by US Patent No. 4568649, European Patent No. 0154734 and Japanese Patent Application No. 84/52452. CyDye or portion thereof is manufactured under licence from Carnegie Mellon University, US Patent Numbers 5,268,486 and 5,569,587 and Patent Application number WO 99/31181.

Cy, CyDye and FARCyte are trademarks of Amersham Biosciences Limited or its subsidiaries.

Amersham Biosciences is a trademark of Amersham plc.

Amersham Biosciences AB SE-751 84 Uppsala Sweden. Amersham Biosciences Inc 800 Centennial Avenue PO Box 1327 Piscataway NJ 08855 USA.

Amersham Biosciences UK Limited Amersham Place Little Chalfont Buckinghamshire England HP7 9NA. Amersham Biosciences Europe GmbH Munzingen

Strasse 9 D-79111 Freiburg Germany. Amersham Biosciences KK Sanken Building 3-25-1 Hyakunincho Shinjuku-ku Tokyo zip 169-0073 Japan

Pefabloc is a trademark of Pentapharm AG, Switzerland; Microbeta is a trademark of Wallac Oy, Finland.

\* To whom all correspondence should be addressed.

©Amersham Biosciences UK Limited, 2001 - All rights reserved.

All goods and services are sold subject to terms and conditions of sale of the company within the Amersham Biosciences

group which supplies them. A copy of these terms and conditions are available on request. Some of these products may

only be available to collaborators and customers within certain of our technology access programmes.

This poster was presented at the 7<sup>th</sup> SBS Annual Conference and Exhibition, Baltimore, MD, USA, 10 - 13 September 2001.