

# NIU Standard Operating Procedure Template



<b>#1</b>	<b>CONTACT INFORMATION: Timothy Hagen, Office: Faraday 350, Phone: (815)-753-1463</b>
<b>Procedure Title</b>	Using the Biotage Initiator Microwave Synthesizer
<b>Procedure Author</b>	Zheng Zhang
<b>Date of Creation/Revision</b>	9/22/2015 (Travis Helgren)
<b>Name of Responsible Person</b>	Timothy Hagen
<b>Location of Procedure</b>	LaTourette Hall (Lab: 334)
<b>Approval Signature</b>	
<b>#2</b>	<b>THIS STANDARD OPERATING PROCEDURE (SOP) IS FOR A:</b>
<input type="checkbox"/> Specific laboratory procedure or experiment <input checked="" type="checkbox"/> Generic laboratory procedure that covers several chemicals <input type="checkbox"/> Generic use of specific chemical or class of chemicals with similar hazards	
<b>#3</b>	<b>PROCESS OR EXPERIMENT DESCRIPTION</b>
The use of the Biotage Initiator Classic for the microwave assisted synthesis of compounds	
<b>Frequency:</b>	<input type="checkbox"/> one time <input type="checkbox"/> daily <input type="checkbox"/> weekly <input type="checkbox"/> monthly <input checked="" type="checkbox"/> other: Based on each experiment
<b>Duration per Expt:</b>	around 5-15 minutes; or _____ hours
<b>#4</b>	<b>SAFETY LITERATURE REVIEW &amp; HAZARD SUMMARY</b>

# NIU Standard Operating Procedure Template



Review MSD sheet for chemicals used in the specific synthesis. Make sure all reactions are kept below 20 bar of pressure.

#5

## STORAGE REQUIREMENTS

Shutdown instrument after use

#6

## STEP-BY-STEP OPERATING PROCEDURE

### Steps to include in your procedure:

1. Don personal protective equipment.

appropriate street clothing (long pants, close-toed shoes)

gloves; indicate type: Nitrile examination gloves

safety goggles  safety glasses  face shield

lab coats

other: \_\_\_\_\_

2. Check the location and accessibility of the safety equipment that serves your lab:

ITEM	STATUS
Laboratory Fume Hood/Glove Box or other Ventilation Control	Location: <u>Various Locations</u>
Eyewash/Safety Shower	Location: <u>Near front lab door/in the hallway</u>

3. Weigh appropriate amounts of each chemical necessary for your reaction into the appropriately sized reaction vial (note: the microwave synthesizer requires the use of specialized vials and lids). Seal the vial using the lid-sealing tool supplied by Biotage. Place the vial in the reaction chamber in the microwave synthesizer.

4. Set your reaction conditions (generally only temperature and time) using the touch screen LED on the front of the instrument. Hit start to run your reaction. Make sure to monitor the pressure inside your reaction vial (the LED screen shows a plot of power output, temperature and pressure). If the pressure exceeds 20 bar, the seals will fail on the microwave synthesizer and an explosion may occur. Therefore, never exceed 18 bar of pressure. Once your reaction runs to completion, the

# NIU Standard Operating Procedure Template



microwave synthesizer will cool your reaction vial to 50 deg C, at which point the sheath will retract and you can remove your vial. Remove the lid using pliers and remove the reaction contents as described in each specific synthetic procedure. Rinse lid with acetone and discard in the trash. Rinse vial with acetone and discard in glass waste container.	
5. Dispose of hazardous solvents, solutions, mixtures, and reaction residues as hazardous waste. See EH&S Hazardous Waste Program <a href="http://www.ehs.niu.edu/ehs/chemical/waste.shtml">http://www.ehs.niu.edu/ehs/chemical/waste.shtml</a>	
6. Clean up work area and lab equipment.	
7. Remove PPE and wash hands.	

<b>#7</b>	<b>WASTE DISPOSAL</b>
-----------	-----------------------

All waste may be placed in the appropriate organic waste container.

<b>#8</b>	<b>TRAINING REQUIREMENTS</b>
-----------	------------------------------

- General Training** *(check all that apply):*
- General Safety & Emergency Preparedness
  - Chemical Safety for Laboratories
  - Radiation Safety
  - Biosafety training
  - Other: Equipment cleaning procedure

<b>Location Where Records Maintained:</b>	
---	--

- Laboratory-specific training** *(check all that apply):*
- Review of SDS for other chemicals involved in process/experiment
  - Review of this SOP
  - Other: \_\_\_\_\_

<b>Location Where Records Maintained:</b>	
---	--

<b>#9</b>	<b>PRIOR APPROVALS</b>
-----------	------------------------

# NIU Standard Operating Procedure Template

---



*Prior approvals are required by the following University Committees:*

*Radiation Safety Committee: Radioactive material,*

<http://www.ehs.niu.edu/ehs/lasersafety/RAM/index.shtml>

*Radiation Safety Committee: X-Ray machines*

<http://www.ehs.niu.edu/ehs/lasersafety/XRay/index.shtml>

*Laser safety: Laser producing equipment Class 3b or above.*

<http://www.ehs.niu.edu/ehs/lasersafety/Laser/index.shtml>

*IACUC: Animal use in research*

[http://www.orc.niu.edu/orc/animal\\_research/index.shtml](http://www.orc.niu.edu/orc/animal_research/index.shtml)

*IBC: Recombinant DNA, potential pathogens, human tissue/body fluids*

<http://www.orc.niu.edu/orc/biosafety/niupolicy.shtml>