

SeqScape Software v2.5

Quick Reference Guide

Applied Biosystems SeqScape Software v2.5 is expressly designed for mutation detection and analysis, SNP discovery and validation, pathogen subtyping, allele identification, and sequence confirmation.

This Quick Reference Guide describes the four workflows performed in these types of studies:



Workflow for Creating a Project



Workflow for Reviewing Data



Workflow for Modifying Settings



Workflow for Exporting and Printing

Workflow for Creating a Project

All analysis in SeqScape software occurs in a project. Before you can create a project you must create a project template that contains:

- A Reference Data Group (RDG)
- Analysis Defaults
- Display Settings

Once a project is created, the project template can be used recurrently.

Note: VariantSEQR customers should proceed directly to Step 4 “Add Sample Files” on Page 7.

Creating a project requires that you follow these steps:

Creating a Project
1. Import the Reference Sequence
2. Define Analysis and Display Settings
3. Create a Project Template
4. Add Sample Files
5. Analyze the Data


Step 1: Import a Reference Sequence

IMPORTANT! The Reference Sequence format can be any of the following:

- .txt, ab1, .seq, .fasta files
- Genbank format files
- Aligned sequences in .fasta (FASTA format)


The workflow described below uses a Genbank file for your reference data. To begin, first download a Genbank file, then create a known variants table.

Note: If you download a Genbank file in the NCBI web interface, select **Display as Genbank**, then **Send to File** to save it as a text file that can be imported.




To use a Genbank file as your reference:

- 1. Open the **SeqScape Manager** from the Tools menu.




- 2. Select the **Reference Data Group** tab.
- 3. Click **New**, then enter a name.



- 4. Select the **ROI (Regions of Interest)** tab.

Note: For information on any of the four panes in the ROI tab (Layer, ROI, Reference Sequence, Reference Segment), click **Info**.

- 5. Click on the **Add Ref. Segment** button.



- 6. Select the Genbank file, then click **Import**.

Note: When the Genbank file imports, information from the feature table populates the Regions of Interest (ROI) table in the SeqScape software and layers are automatically created.

```
FEATURES             Location/Qualifiers
     source            1..792
                        /organism="Homo sapiens"
                        /db_xref="taxon:9606"
                        /map="6p21.3"
     gene              <1..>792
                        /gene="HLA-C"
                        /allele="HLA-Cw*0802"
                        join(<1..270,517..>792)
                        /gene="HLA-C"
                        /product="HLA-C class I antigen"
     mRNA              join(<1..270,517..>792)
                        /gene="HLA-C"
                        /codon_start=3
                        /product="HLA-C class I antigen"
                        /protein_id="AAG53742.1"
                        /db_xref="GI:12483747"
                        /translation="SHSMRYFYTAVSRPRGRGEPRFIAVGYYDDTQVQFDSDAASPRG
EPFAPWVEQEGPEYWDRETQYKRQAQTDVRSLRNLRGYYNQSEAGSHTLQRMYGCDL
GPDGRLLRGYNQFAYDGKDYIALNEDLRSLWTAADKAAQITQRKEWAAREAEQRRAYLE
GTCVEWLRRYLENGKKTQLRA"
     exon              1..270
                        /gene="HLA-C"
                        /number=2
     exon              517..792
                        /gene="HLA-C"
                        /number=3
```

7. Select the **NT Variants** tab (RDG Properties) and select the NT Variants you want to add to the reference sequence, then click **Import** to import a tab-delimited variants file or a multi-aligned sequences (.fasta) file.

Note: When importing an AA variants file, use a tab delimited variants format.

8. Click **OK**.

You have now completed importing the reference sequence.



	A	B	C	D	E	F	G	H
1	Type	ROI	NT position	Reference	Variant	Style	Descr	Used by all ROIs
2	change base	PLAB_target_e_40	40	G	C	Cyan		yes
3	change base	PLAB_target_p_123	123	T	C	Cyan		yes
4	change base	PLAB_target_e_56	56	C	Y	Cyan		yes
5	change base	PLAB_target_p_511	511	T	C	Cyan		yes

RDG Properties								
General ROI NT Variant AA Variants Variant Style								
Type	ROI	Position	Reference	Variant	Style	Description	Used by all R...	
Change Base	PLAB_target_e_40	40	G	C	Cyan		yes	
Change Base	PLAB_target_p_123	123	T	C	Cyan		yes	
Change Base	PLAB_target_e_56	56	C	Y	Cyan		yes	
Change Base	PLAB_target_p_511	511	T	C	Cyan		yes	

Step 2: Define Analysis and Display Settings

To define analysis defaults:

1. Select **Tools > SeqScape Manager**.
2. Select the **Analysis Protocols** tab.
3. Click **New**, then enter a name.
4. Select the **Basecalling** tab, then select the Basecaller file and Dye/Primer file from the drop-down lists.
5. Keep the default settings for Processed Data, Ending Base and Quality Threshold options unless a custom setting is required.
6. Select the **Mixed Bases** tab to define the secondary peak threshold for mixed base identification.



7. Select the **Clear Range** tab.
8. Keep the default settings for Quality Value and Reference Trimming.

Analysis Protocol Editor

General | Basecalling | Mixed Bases | Clear Range | Filter

Clear Range Methods

☐ Use clear range minimum and maximum
First Base >= 20 End Base 550
☐ Bases to trim from 3' end 20

☒ Use quality values
Remove bases from the ends until fewer than 4 bases out of 20 have QVs < 20

☐ Use identification of N calls
Remove bases from the ends until there are fewer than 4 Ns out of 20 bases

☐ Mask M13 universal sequencing primers

☒ Use reference trimming

Multiple clear range methods are applied in order.
Smallest clear range is the result.

OK Cancel

9. Select the **Filter** tab.
Keep the default Filter settings as shown.
10. Click **OK**.
11. Select the **Analysis Defaults** tab in the SeqScape Manager.

Analysis Protocol Editor

General | Basecalling | Mixed Bases | Clear Range | Filter

Filter Settings

Maximum Mixed Bases (%) : 20.0

Maximum Ns (%) : 10.0

Minimum Clear Length (bp) : 50

Minimum Sample Score : 20

☐ Skip filtering if sample-level HIM is detected

12. Click **New**, then enter an Analysis Defaults name.
13. Select the **Sample** tab.

New Analysis Settings

General | Project | Specimen | Sample

Name

Analysis Defaults Name: 3730_AnalysisDefaults

Created: N/A Created By: N/A
Modified: N/A Modified By: N/A
Source: N/A

14. In the drop-down list, select the Analysis Protocol you created.
15. In the **Project** and **Specimen** tabs, keep the default settings.
16. Click **Save**.

New Analysis Settings

General | Project | Specimen | Sample

Settings

Analysis Protocol: 3730_ReSequencingProtocol

☒ Always use this analysis protocol

Generally you will keep the defaults for **Display Settings** and proceed to Step 3, “Creating a Project Template.”

Step 3: Create a Project Template

To create a project template:



1. Select **Tools > SeqScape Manager**.
2. Select the **Project Templates** tab.
3. Click **New**, then enter a Project Templates name.
4. In the drop-down lists, select the Reference Data Group and Analysis Defaults that you previously created in Steps 1 and 2.
5. Keep the default Display Settings.
6. Click **OK**.
7. Close the SeqScape Manager.

You now have created the Project Template.

Step 4: Add Sample Files

To create a project:

1. In SeqScape software, select **File > New Project**.
2. Enter a name for the project.
3. In the Project Template list, select the project template you previously created (in Step 3).
4. Click **New**.



- To add data to your project, select **File > Import Samples**.

IMPORTANT! Groups of sample files from the **same** biological source must be grouped together (into one specimen) in SeqScape software.


- Navigate to your data folder, then group your sample files into specimens.
- For manual specimen creation:
 - Click **New Specimen**.
 - Select files belonging to the specimen, then click **Add Samples**.
- For automatic specimen creation:
 - Define the delimiters by assigning alphanumeric values in the **Start** and **End** fields.
 - Select the data folder, then click **Auto Add**.

Note: It is recommended that you use a sample naming convention such as:

*IndividualID_Gene_Exon_
WellLoc_F/R.ab1*

- Click **OK**.

You now have created a project with samples and you are ready to analyze the data.




Step 5: Analyze the Data

To analyze your data:

1. Click  (Analyze) for analysis to begin.

Note: A progress bar indicates that the software is performing basecalling, trimming, assembly, consensi alignment and reference comparison on your data.

After a successful analysis, your data appear assembled in the Project view.



Workflow for Reviewing Data

Reviewing data require that you follow these steps:


Reviewing Data and Results
1. Review Data per Specimen and Assign Pass/Fail Status
2. Review All Data in Project View
3. View Variant Results in Mutations/AA Variants Reports
4. View Resequencing Results in the Genotyping Report

Step 1: Review Data per Specimen and Assign Pass/Fail Status

To review data by specimen:


1. In the Project Navigator (Project view), select a specimen.

Note: The specimen view shows the sample distribution within the specimen across the reference sequence. Forward and reverse orientations are indicated by arrows representing sequence coverage.




2. In the Project Navigator, open the Analysis QC Report to view failed sample traces (Unassembled).

Note: SeqScape software automatically filters out low quality sample traces based on the settings in the Filter tab of the Analysis Protocol. These samples are grouped in the unassembled node after analysis. Descriptions of the failed samples appear in the Sample Analysis Table (QC Report).



Note: The Analysis QC Report provides a list of sample traces with possible heterozygous indel mutations (HIMs). HIMs are detected when two alleles differ by an insertion or deletion in a sample trace.



Possible Heterozygous Indel Mutations			
Specimen	Sample	Position	Size
Specimen_NA00893	NA00893dup_hCG17690_114_19866955806164.r	307	17
Specimen_NA00893	NA00893dup_hCG17690_114_19866955564566.f	260	17

- Adjust analysis settings to allow failed samples to assemble in the project (see Modifying Settings Workflow for more detail).

Note: You may need to modify your Filter and/or Clear Range analysis settings and then reanalyze data to allow failed samples to pass (Assembled).

Specimen	# Samples	Basecalling	Filter	Assembly	Specimen Score	Total # Variants	Comments
Specimen1	12				40	19	
Specimen2	13				44	15	
Specimen3	14				40	14	
Specimen4	15				43	15	

Specimen 4 has few samples which did not pass through filtering and assembly.

Complete Partial Output No output

Description table gives details of why samples failed filtering and assembly.

Specimen	Sample	Step	Description
Specimen1	Specimen1_PLAB_25_A1.f	Filter	Clear range less than minimum allowed (0-50). Sample score less than minimum allowed (0-25).
Specimen2	Specimen2_PLAB_24_A1.f	Filter	Sample score less than minimum allowed (24-25).
Specimen3	Specimen3_PLAB_25_A1.f	Assembly	Incomplete results presented from previous stage.
Specimen4	Specimen4_PLAB_24_A1.f	Assembly	Incomplete results presented from previous stage.

Filter Settings are set in the Analysis Protocol.

- Select a Segment name, then select the **Assembly** tab to view passed sample traces (Assembled).


Note: The Segment Assembly displays the assembly of sample traces within a specimen for a specific region of the reference sequence. To zoom in/out, select **Cntl +/-**.

- Edit bases with low quality values if necessary.
 - You can customize the bar colors for the quality values (in Display Settings) and flag bases with low quality values in red, for example.




- b. When switching between the Project view and the Assembly view, the base position is bookmarked for easier editing.
- c. Make sure **Discrepancy** is checked in Multiple Tab Jump Settings, if required.

Note: The algorithm automatically corrects for sequencing errors by reviewing data quality in both directions (forward and reverse) to arrive at an accurate consensus.



6. Right-click the specimen name to assign Pass/Fail review status.
7. Repeat Steps 1 to 6 for all specimens in the project.


You are now ready to review your data at the project level.



Step 2: Review all Data in Project View

To review all data:

1. Examine variants between the specimen consensus and the reference sequence by either of the following methods:
 - a. Select the **Character/ Dots** view. In this view, dots appear for all bases as in the reference, and characters appear for variants.
 - b. View data by layers. Select a layer from the **Active Layer** drop-down list.






c. Select the **Tab Jump to Next** field. In this setting, you can select **Known/Unknown Variant**.

2. To add a Known Variant right-click a base, then select **Add Variant**.

3. Click **OK**. The variant is displayed in the Overview pane as a red bar.

4. To view Amino Acid translation, click on the **Amino Acid** view.

You are now ready to review the variant results in reports.



Step 3: View Variant Results in Mutations/AA Variants Reports

To view variant results in reports:

1. Select the desired layer using the **Active Layer** drop-down list, and highlight the project name within the Project Navigator.

Note: To view results at the specimen level, highlight the specimen name.

2. Select **Analysis**, then select **Report Manager**.

3. In the left pane, select the **Mutations Report**.

4. In the right pane, view mutations in the Mutations Table.

5. In the Mutations Table, click on a base change to link the mutation to a specimen consensus in the Project view.

Note: You can review the mutation by looking at its quality value and its position in a region of interest (ROI).



Note: You can sort mutations by quality value, type, and so on by double-clicking any of the column headers. To remove a column, right-click, then deselect the column name.

- Click the triangle icon next to the specimen name to view the sample electropherogram.

Note: Look for mutations at the specimen level or “confirmed” mutations due to bidirectional coverage.

- Open the Mutations report to view the amino acid effect (AA Change/Effect column) for translated regions.
- Click any of the amino acids in the AA Change column to navigate to the AA Variants Report.
- Click an AA Change to link the AA variant to a specimen consensus in the Amino Acid view.

Note: To adjust how the report displays text results, click **Report Settings**, then deselect **Wrap Text**.

To view reports and data simultaneously, select **Window > Tile** in the main menu.

The screenshot shows the 'Mutations' report in a software application. The left pane displays a tree view of the project structure, including 'Project', 'Specimen1', 'Specimen2', and various segments. The main pane shows a table of mutations for Specimen1 and Specimen2. The table has columns for Specimen, Base Change, ROI, Position, Length, Type, QV, Known, and Effect. Below the table, there is a section for 'Click to view electropherogram' showing a sequence alignment and a chromatogram.

Specimen	Base Change	ROI	Position	Length	Type	QV	Known	Effect
Specimen1	400>C	NM_004864_exon_1	400	1	Sub	11	no	missense
Specimen2	400>C	NM_004864_exon_1	400	1	Sub	21	yes	missense
Specimen2	150>S	NM_004864_exon_1	150	1	Sub	12	no	stop
Specimen2	580>V	NM_004864_exon_2	580	1	Sub	24	yes	stop

Specimen	AA Change	Position	Length	Type	Known	NT Change	Description
Specimen1	V14L	14	1	Sub	no	400>C ROI: NM_004864_exon_1	
Specimen1	S53[T,S]	53	1	Sub	no	157T>W ROI: NM_004864_exon_1	
Specimen1	E326Q	326	1	Sub	no	685A>G ROI: NM_004864_exon_2	
Specimen1	D327[D,E]	327	1	Sub	no	689C>S ROI: NM_004864_exon_2	
Specimen2	V14[L,V]	14	1	Sub	no	400>S ROI: NM_004864_exon_1; 400>C ROI: NM_004864_exon_1	
Specimen2	E326Q	326	1	Sub	no	685A>G ROI: NM_004864_exon_2	

The screenshot shows the 'Reports Settings' dialog box. The 'Report Display Settings' tab is selected. The 'Horizontal Scrolling' checkbox is unchecked. The 'Font' is set to Arial. The 'Size' is set to 10. The 'Wrap Text' checkbox is checked.

Step 4: View Resequencing Results in the Genotyping Report

IMPORTANT! The Genotyping Report contains some tables that are filled with the VariantSeqr™ Resequencing System project templates.

To view resequencing results in a report:

- 1. In the Navigation pane, select the project name.

Note: To view results at the specimen level, select the specimen name.

- 2. Select **Analysis**, then select **Report Manager**.
- 3. In the left pane, select **Genotyping Report**.

Specimens in Report

Specimen2

Gene Summary

Gene Symbol

PLAB

Gene Aliases

PDF, MIC1, GDF15, MIC-1, NAG-1, PTGFB, GDF-15

Cytogenetic Band

19p13.1-13.2

Gene Name

prostate differentiation factor

Chromosome

19

LocusLink ID

9518

Transcript Table

Layer

Amplicon ID

Assay Target Start

Assay Target Length

5' Regulatory Length

Exon Length

Intron Length

NM_004864

PLAB_28

12832

300

0

0

0

NM_004864

PLAB_27

12984

436

0

0

0

NM_004864

PLAB_42

13258

355

0

0

0

NM_004864

PLAB_25

13542

454

0

4

0

NM_004864

PLAB_24

13858

421

0

287

0

Resequencing Coverage

Specimen

FOR Cov.

FOR Avg QV

FOR Fwd Cov.

FOR Rev Cov.

5' Reg Cov.

5' Reg Avg QV

5' Reg Fwd Cov.

5' Reg Rev Cov.

Exon Cov.

Exon Avg QV

Exon Fwd Cov.

Exon Rev Cov.

Intron Cov.

Intron Avg QV

Intron Fwd Cov.

Intron Rev Cov.

Specimen2

100.0

46.2

100.0

95.8

0.0

0.0

0.0

0.0

100.0

46.2

100.0

95.8

0.0

0.0

0.0

0.0

Genotype Table

Specimen

T / 157 / NM_004864_exon_1

Specimen2

T (41)

- 4. In the Resequencing Coverage Table, review the resequencing coverage percentage for each specimen.

Note: The resequencing coverage is based on the target regions.

- In the Project view, add genotypes in the Genotyping Report by right-clicking a base, then selecting **Add Genotype**.

Note: To view the genotypes for each specimen, open the Genotyping Report.

- Review the genotypes to compare variants or bases across specimens in a project.

If you completed reviewing your data, and reanalysis is required, go to “Workflow on Modifying Settings” on the next page. If reanalysis is *not* required, you can go directly to “Workflow on Exporting and Printing.”



Workflow for Modifying Settings

You can make modifications to your analyzed project by changing the settings (analysis, reference or display) after reviewing your analysis results.

Modifying settings requires that you follow these steps:

Modifying Analysis Settings
1. Modify Analysis Settings in the Project
2. Modify Reference Settings in the Project
3. Modify Display Settings in the Project
4. Apply a Project Template to an Existing Project

Step 1: Modify Analysis Setting in the Project

To modify analysis settings for *all* samples in a Project or Specimen:



1. Select **Project** or **Specimen Name** in the Project Navigator.
2. Select **Analysis > Sample Manager**.
3. If you want to adjust the Basecaller and Dye Set/Primer file, select the desired settings, then click **Apply**.
4. Click **OK** to close the Sample Manager.

Note: In the project, a red slash through the files that you adjusted indicates the samples are ready to be reanalyzed.

5. Click  (Analyze) to analyze the specimen.


If you need to adjust *more* than the Basecaller and Dye Set/Primer File:

1. Create a new Analysis Protocol in the SeqScape Manager. (See Step 2 on Page 4).



Sample File Name	Specimen	Barcode Name	Basecaller	DyeSet/Primer	Spacing	Peak 1	Start	Stop	Assembled
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	10.77	1497	1535	2692	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.25	1552	1581	7744	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.16	1551	1580	6820	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	10.97	1495	1524	6882	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.25	1541	1570	7885	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.15	1514	1543	7747	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.26	1574	1603	7421	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.22	1557	1586	7584	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.24	1557	1586	7630	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	10.98	1503	1532	7471	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.25	1562	1591	7909	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.19	1538	1567	7981	
Specimen1_FLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3730_POP7_BOTV3.mob	11.04	1509	1538	7774	

2. Select the desired samples in the Sample Manager, then click **Apply** to apply the Analysis Protocol.
3. Click **OK** to close the Sample Manager.



Note: In the project, a red slash through the specimen/ files that you adjusted indicates that the samples are ready to be reanalyzed.

4. Click  (Analyze) to reanalyze the data.

Sample File Name	Specimen	Sample Name	Base Caller	DeSeqFilter	Spacing	Peak 1	Start	Stop	Assembled
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	
Specimen1_PLA_Specimen1	Specimen1	00000019866	KB.bcp	KB_3130_POP7_BDTv3.mob	0.00	0	0	0	


If you are adjusting the Analysis Settings in a Project for just **one** sample:

1. Select the sample in the Project Navigator.
2. Go to **Analysis > Analysis Settings** and make your modifications to the settings.

IMPORTANT! The change applies only to the selected sample and does not modify the Analysis protocol created in the SeqScape Manager.

Note: A red slash through the selected sample indicates the sample is ready to be analyzed.

3. Click on the  icon to reanalyze the data.




Step 2: Modify Reference Settings in the Project

To make changes such as adding variants to an existing reference:

1. Click **rdg** in the project view.
2. Select the **ROI** tab.
3. Modify the reference settings within a project by:
 - Changing the Codon Start Number for amino acid translation
 - Modifying the translation frames
 - Changing the ROI Start numbering in the ROI pane
 - Adding variants and changing variant styles or colors


Note: These modifications are saved only within the *existing* project. They do not change in the original Reference Data Group that you created in the SeqScape Manager.



Step 3: Modify Display Settings in the Project

To modify the Display Settings within a project:


1. In the Project view, click **ds** (Display Settings) to modify viewing options.
2. Select the **Bases** tab.
3. Change the color scale as desired to modify the quality value bars.
4. Select the **Views** tab.



5. Increase the number of bases you want to display in the Project view as indicated in the EP Window field.

6. (Optional) Select **Italicize Reverse Strand**.


Note: The display settings show 70 bases and italicize reverse strands by default.



Step 4: Apply a Project Template to an Existing Project

To make changes to an existing Project Template:


1. Select **Tools > SeqScape Manager**, or simply click **Tooltip**.
2. Modify the data objects within SeqScape Manager and re-select the same data objects within the Project Template to apply changes.
3. Click **Open Folders** to verify your modifications.
4. Click **OK** to close the Project Template editor.



To apply the modified Project Template:

1. Open the project.
2. Select **Apply Project Template** (Analysis menu).
3. Select the appropriate Project Template from the drop-down menu.
4. Click  (Analyze) to reanalyze the data with the modified template.

Note: All previous edits are lost.




Workflow for Exporting and Printing

In this workflow we review optional selections for exporting and printing your analysis results.

Exporting and Printing
1. Select Export Options for an Analyzed Project
2. Set Print Options at Project/Specimen/Segment Levels

To determine which printing and exporting options are available in SeqScape, select a level in the Project Navigator from one of the four choices:

- Project
 - Specimen
 - Reference Segment
 - Sample
1. Select **File > Export, File > Print.**



Step 1: Select Export Options for an Analyzed Project


To export a selected report from a Project:

1. Open the project.
2. Select **Report Manager** (Analysis Menu).
3. Select **Export > Report.**
4. When the Export window opens, select for the format tab-delimited text (Microsoft Excel®), PDF, HTML, or XML.



To export all reports automatically after analysis:

1. Select **Tools > Options**.
2. Select **Export Reports After Analysis**, then in the Format drop-down list, select a format to export.



3. Browse to set a target directory for the file.
4. Click **OK**.


Note: Once a project is reanalyzed, the reports will *automatically* export to the set directory.

To export alignment files:

1. Select the project (Project Navigator), then select **File > Export**. You can export the nucleotides project alignment or the amino acid project alignment.

Note: The reference sequence and the specimen consensus sequences also export.

2. (Optional) Select **File > Export** for the Specimen level/ Reference Segment level to export the consensus sequence or the aligned sample sequence.
3. (Optional) Select **Electropherogram** to export the Sample Sequence.




Step 2: Set Print Options at the Project, Specimen, Segment Levels

Select **All Data** to print the specimen consensus sequences and the reference sequence. Selecting **All Data** does *not* print the Electropherogram data.


Select **Only the Visible Data** to print only visible data. Consequentially, only the *displayed* electropherogram data (specimens) are printed.

Note: The electropherogram data for the full project are *not* printed. Only the data that display in the Project window prints.


From the Specimen view, you can print only the data that are displayed.



From the Segment Assembly view, you can print all data *and* the visible data that are displayed.



By selecting a sample you can print the sample electropherogram.



For more information on Using SeqScape software v2.5, refer to the *SeqScape Software User Manual*.



iScience. To better understand the complex interaction of biological systems, life scientists are developing revolutionary approaches to discovery that unite technology, informatics, and traditional laboratory research. In partnership with our customers, Applied Biosystems provides the innovative products, services, and knowledge resources that make this new, **Integrated Science** possible.

© Copyright 2004, Applied Biosystems. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

Notice to Purchaser: License Disclaimer

Purchase of this software product alone does not imply any license under any process, instrument or other apparatus, system, composition, reagent or kit rights under patent claims owned or otherwise controlled by Applied Biosystems Corporation, either expressly or by estoppel.

SeqScape Software has not undergone specific developmental validation for human identification applications. Human identification laboratories which choose to use SeqScape Software for data analysis should perform their own developmental validation studies.

Applied Biosystems and SeqScape are registered trademarks and AB (Design), Applied and VariantSeq are trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and/or certain other countries. All other trademarks are the sole property of their respective owners.