Miniaturizing Protease Assays using Fluorogenic Substrates and an Optimized Fluorescence Plate Reader

*J. George, R.J. Ward, M.L. Teear, C.G. Norey and D.D. Burns

Amersham Biosciences UK Limited, Amersham Place, Little Chalfont, Buckinghamshire HP7 9NA, U.K. Telephone: + 44 (0)29 2052 6445, Fax: +44 (0)29 2052 6230, e-mail: jeanette.george@eu.amershambiosciences.com

Introduction

Members of the matrix metalloproteinase (MMP) family are responsible for the degradation of components of the extracellular matrix and play an important role in embryonic development, morphogenesis, reproduction and tissue remodelling and repair. Inappropriate MMP activity may contribute to a number of pathologies including arthritis and cardivovascular disease⁽¹⁾.

Fluorogenic peptide substrates provide a particularly convenient enzyme assay method and several laboratories have reported on the design of peptides for monitoring MMP enzyme activity^(2,3). Peptide Mca-RPKPVE[↑](Nva)WRK-Dnp-NH₂ is a well characterised MMP-3 substrate⁽⁴⁾. Previously we labelled this peptide sequence with the donor/quencher pair Cy[™]3/Cy5Q, to give peptide Ac-RPK(Cy3)PVE[↑](Nva)WRK(Cy5Q)-NH₂, and reported the development of a reliable MMP-3 activity assay for use with LEADseeker[™] humogeneous imaging system or FARCyte[™] fluorescence plate reader⁽⁵⁾. Here we show the conversion of this MMP-3 assay from a 384 to a 1536-well format using FARCyte.




Figure 1. Schematic of FRET Protease Assay

Methods

Peptide Substrates: Peptide Ac-RPKPVE(Nva)WRK-NH₂ was synthesised and labelled with active CyDye™ esters using standard methods. Dual labelled peptide was purified by RP-HPLC, purity was assessed by mass spectrometry and material was lyophilised. Prior to use, peptide was reformulated in anhydrous DMSO and stock solutions were stored at -20°C.

Standard Assay: Labelled peptide was incubated with or without recombinant human, catalytic domain MMP-3 (Calbiochem), in 50mM Tris pH7.5, containing 150mM NaCl, 10mM CaCl₂, 10μM ZnCl₂ and 0.05% (w/v) Brij™-35 at 22-25°C. Assays, 384-well format, contained 250nM peptide and 2.5ng/well enzyme in final reaction volumes of 60 μ l. For 1536-well assays, reaction volumes were reduced 10-fold to 6 μ l and peptide and enzyme were added at 125nM and 1.25ng/well respectively. Fluorescence signals were measured on FARCyte using Cy3 filters and dichroic mirrors or on LEADseeker (prototype instrument fitted with a -45°C camera) using 540/570nm excitation and emission wavelengths.

Results

The effect of enzyme concentration on substrate hydrolysis was investigated. The results, presented in figure 2a, demonstrated a linear relationship between fluorescence signal and enzyme concentration over the range 0.02–1.25ng/well. The time course shown in figure 2b was established using 125nM peptide and 1.25ng/well MMP-3. Reactions were complete within 60 minutes at which time we measured >8-fold fluorescence signal increase. No significant signal increase was observed in no enzyme control (NEC) wells.





Figure 2. Enzyme and time dependent hydrolysis of peptide substrate. (a) Peptide (125nM) was incubated with MMP-3 enzyme (0.02 - 10ng/well) as described for 60 minutes. (b) Peptide (125nM) was incubated with 1.25ng/well MMP-3 as described. Data in this and subsequent figures is presented as mean ± 58 M (n = 3). The quality of the miniaturized MMP-3 assay was examined according to the Z'-factor analysis method recently proposed by Zhang and colleagues⁽⁶⁾. The Z'-factor coefficient can be any value \leq 1, and assays with a Z-factor of >0.5 are considered to be reliable and reproducible. The results for the miniaturized MMP-3 assay are presented below. A Z'-factor value of 0.79 was determined.



	Enzyme	Mean Signal	10660
		SD	558.5
	NEC	Mean Signal	1155
		SD	104.1
		S:B	9.2
		Z'	0.79

Figure 3. Z' Factor analysis. Peptide was incubated with (**■**) or without enzyme (**■**) as described for 60 minutes. The experiment has been repeated twice with consistent results. S:B, signal to background ratio.




Figure 4. Inhibitor studies. Peptide was incubated with MMP-3 as described in the presence of 70fM to 3.3 µM inhibitor reformulated in anhydrous DMSO. Fluorescence was measured on FARCyte or LEADseeker as described.

Finally, the effect of the potent, broad range MMP inhibitor N-Hydroxy-1,3-di-(4-methoxybenzenesulfonyl)-5,5-dimethyl-[1,3]-piperazine-2-carboxamide⁽⁷⁾ on peptide hydrolysis was investigated. Figure 4 compares the

results from 384 and 1536-well assays. The IC₅₀ values obtained using the 384-well format were 3nM and 5nM when plates were read on LEADseeker and FARCyte respectively; a value of 11nM was achieved using the 1536-well format FARCyte assay. These values are in good agreement with the reported IC₅₀ for this inhibitor⁽⁷⁾.

Table 1 compares the 384 and 1536-well MMP-3 assays. Miniaturization represented a 10-fold reduction in assay volume, a 20-fold reduction in substrate usage and a 2-fold saving in enzyme costs. Overall, data generated using the different assay formats were comparable.

	384 Assay	1536 Assay	Reduction
Volume	60µl	6μl	10
Peptide	250nM	125nM	20
Enzyme	2.5ng/well	1.25ng/well	2
S:B	11.6:1	8.2:1	-
Z'	0.91	0.79	-
IC 50	5nM	11nM	-

Table 1. Comparison of 384 and 1536 FARCyte assays.

S:B for 384-well format was calculated in optimised assays following digestion of 250nM peptide with 2.5ng/well MMP-3 for 300 minutes. Literature (C₉₆ for inhibitor is 18.4nM⁷⁷)

CONCLUSIONS

- We have successfully miniaturized a 384-well protease assay to a 1536-well format using the FARCyte reader.
- In this case, miniaturization significantly reduced reagent consumption without compromising assay reliability. Miniaturized assay signal-to-background was >8:1 and Z'-factor values of >0.7 were achievable.
- We used MMP-3 as a model protease but the same approach could be applied to other proteases providing suitable peptide sequences are available for labelling with CyDye pairs.

References

- 1. Nagase, H. (1996) Hooper, N.M., ed, pp153-204, Taylor and Francis, London.
- Stack, S.M., and Gray, R.D. (1989) J. Biol. Chem., 264, 4277.
- 3. Netzel-Arnett, S. et al., (1991) Anal. Chem., 195, 86.
- Nagase, H., et al., (1994) J. Biol. Chem., 269, 20952.
 George, J., et al., (2001) Poster, FEBS, Lisbon 2001.
- Beorge, J., et al., (2001) Poster, PEBS, Lisbon 2001.
 Zhang, J-H., et al., (1999) J. Biomol. Screening, 2, 67.
- 7. Pikul, S., *et al.*, (1998) *J. Med. Chem.*, **41**, 3568.

Amersham Biosciences is a trademark of Amersham plc CyDye, Cy, FARCyte and LEADseeker are trademarks of Amersham Biosciences Limited or its subsidiaries. CyDye or

portion thereof le manufactured under levene from Carnegie Mellen University. US Patent Numbers 5,288,466 and 5,560 587 and Patent Application number W0.9931181. CyGO is covered by patent Application number W0.9931181. CyGO is covered by patent Application number 5,288,466 and 5,560 587 and Patent Application number W0.9931181. CyGO is covered by patent Application number N0.9931181. CyGO is covered by patent Application number N0.9931181. CyGO i

Biosciences KK Sanken Building 3:25-1 Hyakunincho Shinjukuku Tokyo 2p 169-0073 Japan Peptide Mca-RPKPVE(Nva)WRK-Dnp-NH₂ is covered by Patent Number 5,770,691. Brij is a trademark of ICI Americas Inc. Mca, 7-methylcoumarin-4-acetic acid; Dnp, 24-dinitrophenyl; ¹, represents enzyme cleavage site All goods and services are sold subject to terms and conditions of sale of the company within the Amersham Bioseinces: group which supplies them. A copy of these terms and conditions are available on request. Some of these products may only be available to collaborators and customers within certain of our technology access programmes. "7 whom all correspondance should be addressed.

