

# Inventory Control Form



## Titanium Cannulated Trochanteric Fixation Nail System

Patient Information:

SYNTHES (USA) To Order: (800) 523-0322  
 SYNTHES (Canada) Ltd. To Order: (800) 668-1119

Date: \_\_\_\_\_

Hospital: \_\_\_\_\_

Surgeon: \_\_\_\_\_

Procedure: \_\_\_\_\_

### Implants

#### Titanium Cannulated Trochanteric Fixation Nails, 10 mm distal diameter, sterile



	RIGHT	LEFT	LENGTH	FEMORAL NECK ANGLE
<input type="checkbox"/>	456.330S	<input type="checkbox"/>	456.331S	300 125°
<input type="checkbox"/>	456.332S	<input type="checkbox"/>	456.333S	320 125°
<input type="checkbox"/>	456.334S	<input type="checkbox"/>	456.335S	340 125°
<input type="checkbox"/>	456.336S	<input type="checkbox"/>	456.337S	360 125°
<input type="checkbox"/>	456.338S	<input type="checkbox"/>	456.339S	380 125°
<input type="checkbox"/>	456.340S	<input type="checkbox"/>	456.341S	400 125°
<input type="checkbox"/>	456.342S	<input type="checkbox"/>	456.343S	420 125°
<input type="checkbox"/>	456.344S	<input type="checkbox"/>	456.345S	440 125°
<input type="checkbox"/>	456.346S	<input type="checkbox"/>	456.347S	460 125°
<input type="checkbox"/>	456.350S	<input type="checkbox"/>	456.351S	300 130°
<input type="checkbox"/>	456.352S	<input type="checkbox"/>	456.353S	320 130°
<input type="checkbox"/>	456.354S	<input type="checkbox"/>	456.355S	340 130°
<input type="checkbox"/>	456.356S	<input type="checkbox"/>	456.357S	360 130°
<input type="checkbox"/>	456.358S	<input type="checkbox"/>	456.359S	380 130°
<input type="checkbox"/>	456.360S	<input type="checkbox"/>	456.361S	400 130°
<input type="checkbox"/>	456.362S	<input type="checkbox"/>	456.363S	420 130°
<input type="checkbox"/>	456.364S	<input type="checkbox"/>	456.365S	440 130°
<input type="checkbox"/>	456.366S	<input type="checkbox"/>	456.367S	460 130°
<input type="checkbox"/>	456.370S	<input type="checkbox"/>	456.371S	300 135°
<input type="checkbox"/>	456.372S	<input type="checkbox"/>	456.373S	320 135°
<input type="checkbox"/>	456.374S	<input type="checkbox"/>	456.375S	340 135°
<input type="checkbox"/>	456.376S	<input type="checkbox"/>	456.377S	360 135°
<input type="checkbox"/>	456.378S	<input type="checkbox"/>	456.379S	380 135°
<input type="checkbox"/>	456.380S	<input type="checkbox"/>	456.381S	400 135°
<input type="checkbox"/>	456.382S	<input type="checkbox"/>	456.383S	420 135°
<input type="checkbox"/>	456.384S	<input type="checkbox"/>	456.385S	440 135°
<input type="checkbox"/>	456.386S	<input type="checkbox"/>	456.387S	460 135°

#### Titanium Cannulated Trochanteric Fixation Nails, 11 mm distal diameter, sterile



	RIGHT	LEFT	LENGTH	FEMORAL NECK ANGLE
<input type="checkbox"/>	456.390S	<input type="checkbox"/>	456.391S	300 125°
<input type="checkbox"/>	456.392S	<input type="checkbox"/>	456.393S	320 125°
<input type="checkbox"/>	456.394S	<input type="checkbox"/>	456.395S	340 125°
<input type="checkbox"/>	456.396S	<input type="checkbox"/>	456.397S	360 125°
<input type="checkbox"/>	456.398S	<input type="checkbox"/>	456.399S	380 125°
<input type="checkbox"/>	456.400S	<input type="checkbox"/>	456.401S	400 125°
<input type="checkbox"/>	456.402S	<input type="checkbox"/>	456.403S	420 125°
<input type="checkbox"/>	456.404S	<input type="checkbox"/>	456.405S	440 125°
<input type="checkbox"/>	456.406S	<input type="checkbox"/>	456.407S	460 125°
<input type="checkbox"/>	456.410S	<input type="checkbox"/>	456.411S	300 130°
<input type="checkbox"/>	456.412S	<input type="checkbox"/>	456.413S	320 130°
<input type="checkbox"/>	456.414S	<input type="checkbox"/>	456.415S	340 130°
<input type="checkbox"/>	456.416S	<input type="checkbox"/>	456.417S	360 130°
<input type="checkbox"/>	456.418S	<input type="checkbox"/>	456.419S	380 130°
<input type="checkbox"/>	456.420S	<input type="checkbox"/>	456.421S	400 130°
<input type="checkbox"/>	456.422S	<input type="checkbox"/>	456.423S	420 130°
<input type="checkbox"/>	456.424S	<input type="checkbox"/>	456.425S	440 130°
<input type="checkbox"/>	456.426S	<input type="checkbox"/>	456.427S	460 130°
<input type="checkbox"/>	456.430S	<input type="checkbox"/>	456.431S	300 135°
<input type="checkbox"/>	456.432S	<input type="checkbox"/>	456.433S	320 135°
<input type="checkbox"/>	456.434S	<input type="checkbox"/>	456.435S	340 135°
<input type="checkbox"/>	456.436S	<input type="checkbox"/>	456.437S	360 135°
<input type="checkbox"/>	456.438S	<input type="checkbox"/>	456.439S	380 135°
<input type="checkbox"/>	456.440S	<input type="checkbox"/>	456.441S	400 135°
<input type="checkbox"/>	456.442S	<input type="checkbox"/>	456.443S	420 135°
<input type="checkbox"/>	456.444S	<input type="checkbox"/>	456.445S	440 135°
<input type="checkbox"/>	456.446S	<input type="checkbox"/>	456.447S	460 135°

### Implants (cont'd)

#### Titanium Cannulated Trochanteric Fixation Nails, 12 mm distal diameter, sterile



RIGHT	LEFT	LENGTH	FEMORAL NECK ANGLE
456.450S	456.451S	300	125°
456.452S	456.453S	320	125°
456.454S	456.455S	340	125°
456.456S	456.457S	360	125°
456.458S	456.459S	380	125°
456.460S	456.461S	400	125°
456.462S	456.463S	420	125°
456.464S	456.465S	440	125°
456.466S	456.467S	460	125°
456.470S	456.471S	300	130°
456.472S	456.473S	320	130°
456.474S	456.475S	340	130°
456.476S	456.477S	360	130°
456.478S	456.479S	380	130°
456.480S	456.481S	400	130°
456.482S	456.483S	420	130°
456.484S	456.485S	440	130°
456.486S	456.487S	460	130°
456.490S	456.491S	300	135°
456.492S	456.493S	320	135°
456.494S	456.495S	340	135°
456.496S	456.497S	360	135°
456.498S	456.499S	380	135°
456.500S	456.501S	400	135°
456.502S	456.503S	420	135°
456.504S	456.505S	440	135°
456.506S	456.507S	460	135°

#### Titanium Cannulated Trochanteric Fixation Nails, 14 mm distal diameter, sterile



RIGHT	LEFT	LENGTH	FEMORAL NECK ANGLE
456.630S	456.631S	300	130°
456.632S	456.633S	320	130°
456.634S	456.635S	340	130°
456.636S	456.637S	360	130°
456.638S	456.639S	380	130°
456.640S	456.641S	400	130°
456.642S	456.643S	420	130°
456.644S	456.645S	440	130°
456.646S	456.647S	460	130°

#### Titanium Cannulated Trochanteric Fixation Nails, 170 mm length, sterile



	DIAMETER	FEMORAL NECK ANGLE
456.314S	10	125°
456.315S	10	130°
456.316S	10	135°
456.317S	11	125°
456.318S	11	130°
456.319S	11	135°
456.321S	12	125°
456.322S	12	130°
456.323S	12	135°

#### Titanium Cannulated Trochanteric Fixation Nails, 235 mm length, sterile



	DIAMETER	FEMORAL NECK ANGLE
456.324S	10	125°
456.325S	10	130°
456.326S	10	135°
456.327S	11	125°
456.328S	11	130°
456.329S	11	135°
456.510S	12	125°
456.511S	12	130°
456.512S	12	135°

#### 11.0 mm Titanium Helical Blades (gold)



LENGTH	LENGTH
456.301° 80 mm	456.300S 75 mm
456.302° 85 mm	456.310S 125 mm
456.303° 90 mm	456.650S 130 mm
456.304° 95 mm	
456.305° 100 mm	
456.306° 105 mm	
456.307° 110 mm	
456.308° 115 mm	
456.309° 120 mm	

#### 11.0 mm Titanium Screws, sterile (gold)



LENGTH	LENGTH
04.032.070S 70 mm	04.032.105S 105 mm
04.032.075S 75 mm	04.032.110S 110 mm
04.032.080S 80 mm	04.032.115S 115 mm
04.032.085S 85 mm	04.032.120S 120 mm
04.032.090S 90 mm	04.032.125S 125 mm
04.032.095S 95 mm	04.032.130S 130 mm
04.032.100S 100 mm	

# Inventory Control Form

## Titanium Cannulated Trochanteric Fixation Nail System



### Implants (cont'd)

#### 5.0 mm Titanium Locking Screws<sup>◇</sup> (green) with 3.5 mm hexagonal drive



	LENGTH		LENGTH		
<input type="checkbox"/>	458.926	26 mm	<input type="checkbox"/>	458.954	54 mm
<input type="checkbox"/>	458.928	28 mm	<input type="checkbox"/>	458.956	56 mm
<input type="checkbox"/>	458.930	30 mm	<input type="checkbox"/>	458.958	58 mm
<input type="checkbox"/>	458.932	32 mm	<input type="checkbox"/>	458.960	60 mm
<input type="checkbox"/>	458.934	34 mm	<input type="checkbox"/>	458.964	64 mm
<input type="checkbox"/>	458.936	36 mm	<input type="checkbox"/>	458.968	68 mm
<input type="checkbox"/>	458.938	38 mm	<input type="checkbox"/>	458.972	72 mm
<input type="checkbox"/>	458.940	40 mm	<input type="checkbox"/>	458.976	76 mm
<input type="checkbox"/>	458.942	42 mm	<input type="checkbox"/>	458.980	80 mm
<input type="checkbox"/>	458.944	44 mm	<input type="checkbox"/>	458.985	85 mm
<input type="checkbox"/>	458.946	46 mm	<input type="checkbox"/>	458.990	90 mm
<input type="checkbox"/>	458.948	48 mm	<input type="checkbox"/>	458.995	95 mm
<input type="checkbox"/>	458.950	50 mm	<input type="checkbox"/>	458.999	100 mm
<input type="checkbox"/>	458.952	52 mm			

#### 4.9 mm Titanium Locking Bolts<sup>◇</sup> (green) with 3.5 mm hexagonal drive



	LENGTH		LENGTH		
<input type="checkbox"/>	459.26	26 mm	<input type="checkbox"/>	459.54	54 mm
<input type="checkbox"/>	459.28	28 mm	<input type="checkbox"/>	459.56	56 mm
<input type="checkbox"/>	459.30	30 mm	<input type="checkbox"/>	459.58	58 mm
<input type="checkbox"/>	459.32	32 mm	<input type="checkbox"/>	459.60	60 mm
<input type="checkbox"/>	459.34	34 mm	<input type="checkbox"/>	459.64	64 mm
<input type="checkbox"/>	459.36	36 mm	<input type="checkbox"/>	459.68	68 mm
<input type="checkbox"/>	459.38	38 mm	<input type="checkbox"/>	459.72	72 mm
<input type="checkbox"/>	459.40	40 mm	<input type="checkbox"/>	459.76	76 mm
<input type="checkbox"/>	459.42	42 mm	<input type="checkbox"/>	459.80	80 mm
<input type="checkbox"/>	459.44	44 mm	<input type="checkbox"/>	459.85	85 mm
<input type="checkbox"/>	459.46	46 mm	<input type="checkbox"/>	459.90	90 mm
<input type="checkbox"/>	459.48	48 mm	<input type="checkbox"/>	459.95	95 mm
<input type="checkbox"/>	459.50	50 mm	<input type="checkbox"/>	459.100	100 mm
<input type="checkbox"/>	459.52	52 mm			

#### Titanium End Caps<sup>◇</sup> (green)



	EXTENSION	
<input type="checkbox"/>	456.311	0 mm
<input type="checkbox"/>	456.312	5 mm
<input type="checkbox"/>	456.313	10 mm

#### Set Screw for TFN, sterile



<input type="checkbox"/>	04.032.000S	125°
<input type="checkbox"/>	04.032.001S	130°
<input type="checkbox"/>	04.032.002S	135°

### Instruments

#### Guide Wire

357.399 3.2 mm Guide Wire, 400 mm

#### Cannulated Drill Bits

357.394 17.0 mm diameter, large quick coupling, 300 mm

357.403 6.0 mm/10.0 mm diameter, stepped, large quick coupling, 435 mm

357.404 11.0 mm diameter, tapered, large quick coupling, 280 mm

#### 4.0 mm Three-Fluted Drill Bits, quick coupling

315.40<sup>◇</sup> 195 mm  
 357.407 260 mm, 65 mm calibration  
 511.417 150 mm, brad point

#### Reaming and Guide Rods

351.71<sup>◇</sup> 3.0 mm Reaming Rod, 950 mm, with offset ball tip  
 351.76<sup>◇</sup> 3.0 mm Reaming Rod, 950 mm, with straight ball tip

351.704S 2.5 mm Reaming Rod with ball tip and extension, 1150 mm, sterile  
 351.706S 2.5 mm Reaming Rod with ball tip, 950 mm, sterile  
 351.707S 2.5 mm Reaming Rod with ball tip and extension, 950 mm, sterile  
 351.708S 2.5 mm Reaming Rod with ball tip, 1150 mm, sterile

355.041<sup>◇</sup> 3.0 mm Guide Rod, 950 mm, with smooth tip

<sup>◇</sup> Available nonsterile or sterile-packed. Add "S" to product number for sterile product.