Invacare® Top End® T3™ Tennis Elite W heelchair

Price List and Order Form Effective October 1, 2001 Revised May 1, 2002



For ease of ordering, contact Top End Toll Free at: 800-532-8677or 727-522-8677 Fax: 727-522-1007

To order Literature: 800 -828 - 6282

www.invacare.com



			Top End Job Number
Account # _	P.O. #	Phone #	Fax #
Address		City	State Zip Code
Customer N	ame		Phone No
Address		City	State Zip Code
Ship order to: ☐ Provider ☐ Customer Order will ship to Provider unless otherwise noted.			
Chair Type ☐ T3ELT	Invacare Top End Elite 3 wheel custom aluimum f	rame2 305 00	Invacare Top End T-3 Tennis Elite
☐ T4ELT Seat W idth ☐ T3EL23 Seat Depth	Invacare Top End Elite 4 wheel custom aluimun fr	ame2,495.00 "Standard125.00 StandardStandardStandard	
	(13" - 25 - 520 or 25 - 540 wheel size only Adjustable Rear Seat Height* ☐ 13.5" - 17" (*Must order T3EL1 Adjustable Axle Camber Tube ght	125.00 with this option)	
or	☐ Fixed ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ Fixed w/Adjustable Height Range	Standard Standard 20"	Fixed - * Anti-Tip is Optional
w/Fixed Car	I ☐ Fixed (choose one) ☐ 2" ☐ 3" ☐ 4" ☐ 5" ☐ 6 nber (degrees choose one) ☐ 12 ☐ 15 ☐ 18 ☐ 20	5" □ 7"Standard Standard	Back Height Seat Height Rear Aske Position
	Adjustable Axle Camber Tube 2000 (degrees choose one) 12		4. W heel Base 5. Seat Height Front weight limit - 250 lb.
	Add. Camber Tube 2000 (degrees - choose 1 angle) 0 0 3 0 6 0 9 0 12 0 15 18 0 20		Performance W heels (Double wall rim/high flange hub/high pressure tires)
Footrests Mu	ulti-Adjustable System w/ Strap Toe-Clip Cage	Standard	Rear W heel SizeStandard □ 22" (25 - 520) □ 24" (25 - 540) □ 25" (25 - 559) □ 26" (25 - 590) □ 700C □ T3EL31 Spinergy® Spox W heels550.00
Frame Color Standard ColorsStandard			☐ 24" (25 - 540) ☐ 25" (25 - 559) ☐ 26" (25 - 590) ☐ Black ☐ Yellow
ū	☐ Red ☐ SunnyYellow ☐ Cobalt Blue☐ Purple☐ ☐ True Blue☐ ☐ White☐ Anthracite☐ Denim Blue☐ Jewel Green		☐ T3EL3 O mit Build for150.00
□ P 1	TwoTone (choose from standard colors only)	75.00	(Not available on Spingery wheels) T3EL17 KIK® "no flat" Tires75.00
□ P7	Front Back Rugged Color	150.00	☐ Gray ☐ Black ☐ T3EL9 Continental Tubular Tires125.00
□ P19	☐ Silver Vein Candy Color ☐ Electric Red ☐ Electric Blue	200.00	(25 - 540 only) ☐ T3EL8 Spoke Guards75.00
□ P17	☐ Electric Red ☐ Electric Blue ☐ Electric Purple ☐ Electric Teal Sparkle Color (over Black Base Coat only)	200.00	Anodized PushrimsStandard □ Standard □ Close Mount
□ P23	Black w/Twilight Afflair Color	200.00	☐ T3EL29 Coated Pushrims100.00 ☐ T3EL35 DualGrip™ Pushrims75.00
	☐ Midnight Blue		Front W heel SizeStandard
Special Instr	uctions		Options ☐ T3EL5 Cushion 2"

Invacare® Top End® T-3™ Tennis Elite W heelchair Measuring Guidelines

Seat Width - Measure frame outside to outside in 1" increments. Available from 12" to 19" wide.

Tapered Seat - Measure frame width desired at back post to 10" deep and then front at edge of seat upholstery. Maximum difference is 4".

Seat Depth - Measured by the length of upholstery from the front of the seat sling to the back post. Standard is 15" - but is available from 12" to 16".

Seat Rail Extensions - Fixed. Standard length is 3" and is measured from the edge of the seat upholstery.

Seat Height Front / Rear - Floor to top of tube in the front and rear respectively, front range is 17" - 21", rear 13" - 17". Most daily used chairs have a 1" - 3" difference. Tennis players use 4" - 6" difference between front and rear. High level injuries prefer more drop in the seat and lower injuries less. The drop adds balance, security in the seat and better positioning for high injuries. This is very evident in all contact sports such as quad rugby and basketball.

Adjustable Rear Seat Height - An optional 13.5" - 17" adjustable rear seat is available for those who are unsure.

Backrest Height - Ranges from 8" to 16". The back is fixed and it can be made to the exact height or any adjustable range of 11" - 14", 14" - 18" or 16" - 20".

W heelbase - Standard for Tennis is 24". 24" is recommended as a minimum for tennis players who sit above 16" in the rear and for those needing stability, extra leg room and for daily use. The wheelbase is measured from the center of the axle to the center of the front fork. To measure the wheelbase on an adjustable axle position T-3, the axle should be 3" from the front of the back post.

Axle Position - Available as fixed w/fixed camber at 12°, 15°, or 20° for elite players. Only one camber is available w/fixed axle choice or pick. The fixed axle is measured from the front of the back post to center of axle. Minor adjustments can be made w/adjustable back upholstery and foot placement.

Adjustable Axle Position Camber Tube - The best choice for most players, especially beginners or intermediate players who may want to change their axle position as their skills improve. Camber is available in 12°, 15° or 20°. An additional tube is available in camber from 0° - 20°.

Footrest - The footplate system is adjustable both forward and backwards. The dual footplates are multi-adjustable and also can be angled for full foot contact and comfort. Toe strap is standard. Toe-clip cage is available as an option.

Frame W heel Clearance - 1" for fixed axle and adjustable w/adjustable axle camber tube.

Rear W heels - Standard size is 24" (25 - 540). For short players or children, the 22" (25 - 520) wheel is a good choice. The 25" (25 - 559) wheel is good for taller players who want to "push more wheel". Double wall rim and high flange hubs are standard which are great for sports enthusiasts. High pressure tires are standard but threaded tires and KIK tires are also available.

Front Caster - For court sports, the 3" in-line wheel is preferred. 5" is available for daily use.

Options - Sport and everyday options are available.

Color Selection - See Color Selection Guide for details.



Invacare Corporation

www.invacare.com

INVACARE TOP END 4501 63rd Circle North Pinellas Park, FL USA 33781-5914 800-532-8677 Invacare and Top End are registered trademarks and T-3 and Yes, you can. are trademarks of Invacare Corporation. Spinergy is a registered trademark of Spinergy, Inc. KIK is a registered trademark of KIK Technology, Inc. DualGrip is a registered trademark of Sun Metal Products, Inc. Specifications and prices are subject to change without notification.

© 2002 Invacare Corporation

Form No. 96-264 Rev. 05/02