




# FarrowWrap® Thigh High OTS and Custom Order Form

Business: \_\_\_\_\_ Patient: \_\_\_\_\_ PO: \_\_\_\_\_  
 Person Measuring: \_\_\_\_\_ Phone/Fax: \_\_\_\_\_ / \_\_\_\_\_ Date: \_\_\_\_\_



OTS Components		Thighpiece (Only comes with wrap-around kneepiece)									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(G) Groin		60–70 cm		65–75 cm		70–80 cm		75–90 cm		85–100 cm	
(E) Above Knee		47–52 cm		53–59 cm		60–66 cm		67–73 cm		74–80 cm	
(E-G) Short	FW_-O-TS	18 cm		18 cm		18 cm		18 cm		18 cm	
(E-G) Regular	FW_-O-TR	23 cm		23 cm		23 cm		23 cm		23 cm	
(E-G) Tall	FW_-O-TT	29 cm		29 cm		29 cm		29 cm		29 cm	
		Legpiece									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(C) Mid-Calf		36–43 cm		42–50 cm		48–58 cm		53–63 cm		58–68 cm	
(B) Ankle		21–25 cm		25–30 cm		30–36 cm		36–42 cm		42–50 cm	
(A-D) Regular	FW_-O-LR	35 cm		37 cm		39 cm		41 cm		41 cm	
(B-D)		28 cm		30 cm		31 cm		33 cm		33 cm	
(A-D) Tall	FW_-O-LT	40 cm		42 cm		44 cm		46 cm		46 cm	
(B-D)		33 cm		35 cm		36 cm		38 cm		38 cm	
		Footpiece									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(A') Mid-Foot		22–24 cm		25–27 cm		28–30 cm		31–34 cm		35–40 cm	
(X) Regular	FW_-O-FR	16.5 cm		18 cm		19.5 cm		21 cm		23 cm	
(X) Long	FW_-O-FL	18.5 cm		20 cm		21.5 cm		23 cm		25 cm	

Options	
FABRIC	<input checked="" type="checkbox"/>
Classic	<input type="checkbox"/>
LITE	<input type="checkbox"/>
STRONG	<input type="checkbox"/>

Use Nonstandard Liner Option (see liner form)	<input type="checkbox"/>
-----------------------------------------------	--------------------------

Custom Component	Thighpiece	<input type="checkbox"/> Wrap-around <input type="checkbox"/> Standard	Legpiece	FW_-C-L	Footpiece	FW_-C-F	Extend footpiece length to base of toes
Quantity	Right: <input type="text"/>	Left: <input type="text"/>	Right: <input type="text"/>	Left: <input type="text"/>	Right: <input type="text"/>	Left: <input type="text"/>	L   R   BL: <input type="text"/>



## Directions for Thigh High Garments

Only circumferences at points **A**, **B**, **C**, **E**, & **G** and length measurements **X**, **A-D**, & **E-G** are necessary to determine if the patient can fit into an **OTS (Off-The-Shelf or ready-made)** garment. Please note that the **E-G** length measurement, used to determine the length for an OTS thighpiece, can be measured directly (for an OTS garment) or determined by subtracting (A-E) from (A-G), both of which are necessary for a custom thigh high garment. Furthermore, all length measurements are measured along the posterior contour of the leg except A-G, which is measured from the posterior heel/floor intersection to the medial groin, and following contours around any lobes. If any lobes are present, drawings and/or pictures (with the patient's permission) will help. Follow these instructions on how to properly measure.

### FarrowWrap® Thighpiece

1. Obtain the following 3 measurements:
  - d. Upper thigh circumference (at **G**), just below the groin.
  - e. Distal thigh circumference (at **E**), right above the knee.
  - f. Length of thigh from just above the knee joint to the groin. Measure medially. **E-G**.
2. Now, match the OTS upper and distal thigh circumferences to a respective size. If the upper thigh AND distal thigh circumferences do not both match a single size, or if the patient has a large lobule, then a custom thighpiece is necessary along with all measurements.
3. If the circumferences do match an OTS size, then compare **E-G** to the OTS short, regular, and tall length measurements. If **E-G** falls within 2 centimeters of one these lengths, then in the space next to the appropriate OTS length enter the number of thighpieces you would like in that given size. If the measured length is not within 2 centimeters of any of these lengths, a custom thighpiece is highly recommended along with all associated measurements. Enter the number of thighpieces desired.
4. All OTS thighpieces come with wrap-around kneepieces, which is a LITE short-stretch wrap attached to the distal thighpiece. Custom thighpieces come with a wrap-around kneepiece or a standard kneepiece. See our catalog for descriptions. If choosing a custom thighpiece, check the custom thigh sku with the kneepiece you desire, FW\_-C-TWA for wrap-around and FW\_-C-TST for standard.
5. We recommend biker's shorts or GarmentGrip™ to help hold thigh pieces in place on problematic patients.

### FarrowWrap® Legpiece and Footpiece

1. See FarrowWrap® A-D (Knee High) Order Form Instructions

### Note

The measurements used by this form allow for at least a 20% additional reduction in circumference (additional edema reduction), but can only accommodate about a 5% increase. In general, Classic will accommodate greater reductions and LITE greater increases. Patients may fit into an OTS garment even though this form may indicate otherwise, but we cannot guarantee it.

### Liner Options

We recommend at least 2 pairs of liners per extremity. TG® Soft Liners come in 3 circumferences, S <40 cm, M 40-70 cm, & L 70-125 cm.


### Other Options

Hook Stays can be used to help prevent FarrowWrap® legpieces from sliding on highly triangular legs and can be easily attached to the garment. Stays are available in small (16 cm length), medium (21 cm length), and large (27 cm length).



# FarrowWrap® Liners Order Form

Business: \_\_\_\_\_ Patient: \_\_\_\_\_ PO: \_\_\_\_\_  
 Person Measuring: \_\_\_\_\_ Phone/Fax: \_\_\_\_\_ / \_\_\_\_\_ Date: \_\_\_\_\_



Farrow Silver Liners		Sizes			
		Small (Flesh Color)		Large (Silver Color)	
Calf (C) or Arm Circumference (F <sup>1</sup> )		15 – 58 cm		52 – 70 cm	
Silver AD Liner	FLS_-AD_	Qty		Qty	
Silver Arm Liner	FLSF-A2	Qty			
Silver Arm Hand Liner	FLSF-AH2	Qty			


**TG® Soft** is used on patients who are allergic or adversely react to silver, do not fit into a silver liner, who need a little more cushioning than our silver liners, or who desire an open toe.

TG® Soft Liners		Sizes					
		Small		Medium		Large	
Extremity Circumference		< 40 cm		40 – 70 cm		70 – 120 cm	
Arm (F <sup>1</sup> )	FLTG-ARM	Qty		Qty		Qty	
AD Knee High (C)	FLTG-AD	Qty		Qty		Qty	
Are the AD liners a SUBSTITUTE to our standard liners? (Only if you circle "YES" here will you not be charged.)						YES / NO	
AG Thigh High (G)	FLTG-AG	Qty		Qty		Qty	



# FarrowFoam™ Liners Order Form

Business: \_\_\_\_\_ Patient: \_\_\_\_\_ PO: \_\_\_\_\_  
 Person Measuring: \_\_\_\_\_ Phone/Fax: \_\_\_\_\_ / \_\_\_\_\_ Date: \_\_\_\_\_



**FarrowFoam™ Liners** are used under a Farrow brand compression garment. When ordering a FarrowFoam™ Liner you **MUST** add 10 cm to all leg circumferences and 5 cm to foot circumferences when choosing the correct-sized Farrow leg or foot garment to go over it. On the order form, please indicate the actual measurements and then the measurements needed to fit over the FarrowFoam™ liner. XL and XXL sizes are subject to price increases.

FarrowFoam™ Smooth Foam AD Liner		Sizes											
		XSmall		Small		Medium		Large		XLarge*		XXLarge*	
(D) Below Knee		37 – 43 cm		42 – 48 cm		45 – 51 cm		50 – 56 cm		55 – 61 cm		60 – 66 cm	
(C) Mid calf		32 – 38 cm		37 – 43 cm		42 – 48 cm		47 – 53 cm		52 – 58 cm		57 – 63 cm	
(B) Least Ankle		18 – 23 cm		21 – 26 cm		23 – 28 cm		26 – 31 cm		28 – 33 cm		31 – 36 cm	
(AD) Length		38 cm		39 cm		41 cm		42 cm		43 cm		45 cm	
Side Closure AD Liner	FLFM-AD	Qty		Qty		Qty		Qty		Qty		Qty	

© 2013 FarrowMed LLC. All Rights Reserved