

## Features and Fixes for SRS EHR v9.6.029 (Teal)

TTID	Type	Functional Area	Summary
-	Enhanced Feature	Charge Passage	Charge Passage will now include HCPCS codes, modifiers and units.
-	Enhanced Feature	CPOE	Default Future Order Attributes in Order Entry, Favorite Completed Order List and Orders Admin PowerTab to create default future orders.
-	Enhanced Feature	CQM Report	Includes CQM updates for 2015 PQRS and MU reporting.
-	Enhanced Feature	Reports	Provides the ability to export the Referral Aging Report to a CSV file.
10979	Unexpected Behavior	CPOE	Manual lab results entered into SRS now follow the Provider Routing rules and properly display the message.
10411	Unexpected Behavior	EDS	Users will now be able to access the EDS server even if the connection to SRSAPPSVR is lost.
11359	Unexpected Behavior	Labs	Lab results will no longer display the entire requisition in red if only one value is abnormal.
11110	Unexpected Behavior	Reports	Imaging Results measure on the Meaningful Use 2 Report Card will now increment the numerator when manually matching a file to a radiology order via Scan Place.
11503	Unexpected Behavior	Reports	Imaging Results measure on the Meaningful Use 2 Report Card will now increase the numerator for results received via a bi-directional PACS interface.
11283	Unexpected Behavior	Reports	Optimized performance for the Referral Aging report. It no longer freezes or crashes the SRS EHR.
11288	Unexpected Behavior	SRS Portal	A Communicator message will no longer fail on the third response within the same message thread.

## Features and Fixes for SRS EHR v9.5.520 Cloud

TTID	Type	Functional Area	Summary
11412	Unexpected Behavior	eRx / Med History	Medication History requests will now process successfully if a patient has a hyphen present in their 9-digit zip code.
11197	Unexpected Behavior	eRx / Med History	Medication History requests will no longer fail when a patient has either the ARGS or ARGSDIR pharmacy benefit management company (PBM).
11486	Unexpected Behavior	eRx / Med History	Medication History requests will no longer display an error message for patients with an inactive PBM and successfully process the request.
11492	Unexpected Behavior	eRx / Med History	Medication History requests will now process successfully if the patient's PBM data has no start or end date.

## Features and Fixes for SRS EHR v9.5.519 (Jupiter)

TTID	Type	Functional Area	Summary
-	New Feature	Messaging	Emdeon lab results messages now give you the option of having one message per requisition or one message per test ordered. This can be configured for each individual lab for which the practice is receiving results.
-	Enhanced Feature	Charge Passage	When submitting a billing file, we have prioritized E&M codes to be ordered first in the PM system, followed by CPT codes.
11290	Unexpected Behavior	General	Introduced a Workstation Property within the Admin Console to address a race condition when switching too quickly within specific modules that sporadically caused gray lines throughout the application and required the user to re-login to SRS EHR.
11332	Unexpected Behavior	Messaging	Unmatched results previously had messages closed when matched to a patient/order with no forwarded message. Now, when unmatched results are matched to a patient/order, a message is sent to the appropriate person/pool as if the results had come in as matched according to the routing rules defined for that interface.
11289	Unexpected Behavior	Printing	The second and subsequent print jobs now print without errors.

## Features and Fixes for SRS EHR v9.5.025 (Marsala)

TTID	Type	Functional Area	Summary
-	New Feature	Charge Passage	SRS Integration with PM systems allowing for Outbound HL7; including the following message types: SIU, ADT, and DFT.
-	New Feature	CPOE	Certification for LabCorp interface via Emdeon.
10982	Enhanced Feature	Advanced Reporting	Added E-mail and Phone Number filters in the Advanced Reporting Module within the Demographics tab.
10346	Unexpected Behavior	Advanced Reporting	The Advanced Reporting Module Reminders filters are now working as expected when running a patient-reminder report.
10440	Unexpected Behavior	Communicator	Users can now select the Ask Me Later option after clicking on What Is This on a pending chart without causing a login failure.
10848	Unexpected Behavior	Communicator	Using the Chrome browser no longer results in a loop when using SRS Communicator.
10933	Unexpected Behavior	Communicator	Logging in after previous log off from SRS Communicator while using Chrome or Firefox is now permitted.
11275	Unexpected Behavior	Integration	Optimized Integration processing bandwidth and corrected error handling for HL7 files missing insurance fields.
10891	Unexpected Behavior	MU Reports	Users can now export the Drilldown Report after running the MU Report Card.
10892	Unexpected Behavior	MU Reports	Imaging Results excludes "Performed, Results Not Required" and "Performed with Adverse Event" from the numerator in the MU Reports.

## Features and Fixes for SRS EHR v9.5.025 (Marsala) (continued)

TTID	Type	Functional Area	Summary
10932	Unexpected Behavior	MU Reports	The Encounter Filter Dates within the CQM Report are now working as expected.
10934	Unexpected Behavior	MU Reports	Stratification calculations for CMS 155 have been corrected.
10961	Unexpected Behavior	Performance	Optimized SQL Server query that minimizes CPU utilization on the SQL Server, commonly known as the Primary SRS Server.
10862	Unexpected Behavior	Printing	User permissions no longer randomly change when users print from SRS EHR.
10785	Unexpected Behavior	Rx	Nightly automatic retrieval of the patients' Medication History for next day encounters now processes after the 10.6 Surescripts Pipeline update. Additionally, users no longer require a unique NPI in order to request patient Medication History.
10981	Unexpected Behavior	Rx	Corrected missing NDCs from the Rx Drug Database that affected clients on the 10.6 Surescripts Pipeline.

## Features and Fixes for SRS EHR v9.4.019 (Olive)

TTID	Type	Functional Area	Summary
-	New Feature	eRx / EPCS	Electronic Prescribing of Controlled Substances (EPCS) for NY State only.
10353	Unexpected Behavior	Aged Message Report	Users can sign off on Secure messages sent from a Communicator portal account to SRS EHR.
10738	Unexpected Behavior	AMC Report	Not Performed Statuses are no longer included in the numerator for the AMC Image Results measure.
10695	Unexpected Behavior	Communicator	Patients using the SRS Communicator portal can now respond to a non-provider-initiated secure message.
10696	Unexpected Behavior	Communicator	The Communicator temporary password now prints out the correct encounter location or the practice default contained in the SRS EHR Admin Module.
10724	Unexpected Behavior	CPOE	Special instructions are no longer truncating after 250 characters.
10892	Unexpected Behavior	CPOE	Radiology Orders without valid results are correctly calculated and not included in the Imaging Results Menu Measure.
10776	Unexpected Behavior	Direct	Unmatched Direct Messages no longer results in an error message when attempting to match to a patient chart.
10723	Unexpected Behavior	Document Management	Printing no longer randomly overwrites/removes a user's SRS printing permissions.
10725	Unexpected Behavior	Dx App	Corrected the ICD10 drilldown options font from being too light in the Diagnosis Web Application.
10721	Unexpected Behavior	EDS	Users can now proceed with login to an Emergency Document Server if the Primary SRS Server is unavailable.
10688	Unexpected Behavior	Emdeon Labs	Unmatched Emdeon lab results do not receive an error message when attempting to match to a patient chart.

## Features and Fixes for SRS EHR v9.4.019 (Olive) (continued)

TTID	Type	Functional Area	Summary
10722	Unexpected Behavior	Immunizations	Updated VIS dates for various immunizations.
10687	Unexpected Behavior	Labs	Viewing lab results via message center no longer results in error message under specific data conditions.
10707	Unexpected Behavior	Labs	Lab Order Matching no longer results in error message under specific data conditions.
10769	Enhanced Feature	Labs	Service Facility and Report Status are now required fields when manually entering a lab result.

## Features and Fixes for SRS EHR v9.3.025 (Gold)

TTID	Type	Functional Area	Summary
10488	Unexpected Behavior	AMC / Labs	The denominator and numerator for the MU CPOE measure now increments after adding a result via a manual lab entry.
10326	Unexpected Behavior	Document Management	Corrected issue with users being unable to print a specified page range from a multi-page document.
9983	Unexpected Behavior	Forms	Forms can now be edited when patients have results in the Lab Management Module without an ordering provider.
10561	Unexpected Behavior	Radiology	SRS Bi-Directional Radiology for Non-MedStrat interfaces (i.e., Viztek, Merge, Sectra, etc.) now displays the Image Result within SRS EHR.
10562	Unexpected Behavior	Radiology	SRS Bi-Directional Radiology for Non-MedStrat interfaces (i.e., Viztek, Merge, Sectra, etc.) now automatically matches patient's Image Results.
10567	Unexpected Behavior	Source of Payment	Corrected issue with saving patient's demographics as the Source of Payment for Medicare SOP Code 1.
8205	Unexpected Behavior	Source of Payment	Within the patient's demographics the user can now click Save without having a Source of Payment selected.
10595	Unexpected Behavior	CPOE	The CPOE / Orders and Results App now displays the History screen without an error message.
10579	Unexpected Behavior	CPOE Labs	The Emdeon Lab Results interface displays without a 404 Error.
10596	Unexpected Behavior	CPOE Labs	Users can now save manually entered Lab Results.
9423	Unexpected Behavior	CPOE Labs	Lab Reports results now display the correct page number in the footer of the page and include the corresponding Performing Lab Information at the bottom of the result.

## Features and Fixes for SRS EHR v9.3.023 Cloud

TTID	Type	Functional Area	Summary
-	New Feature	CQM / QRDA	Implemented the QRDA Download PowerTab, which provides a link to open/save the file required for electronic reporting.
-	Enhanced Feature	CQM / QRDA	Enhanced and updated the Clinical Quality Measures (CQM) Report interface to an asynchronous report delivery system. This provides the ability for the user to complete other tasks while the report is being processed. An e-mail notification is sent when the final report is ready.
-	Enhanced Feature	CQM / QRDA	Enhanced user interface for the CQM / CDS Admin PowerTab facilitating the ability to electronically report on PQRS.

## Features and Fixes for SRS EHR v9.2.083 (Silver)

TTID	Type	Functional Area	Summary
-	New Feature	CQM / QRDA	Implemented QRDA functionality to facilitate the ability to electronically report on PQRS.
10230	Unexpected Behavior	CCDA / Smoking	Smoking Status within the CCDA was utilizing an incorrect encounter start time, causing CQM calculations to be off.
10460	Unexpected Behavior	CCDA	If a patient encounter did not have a start and stop time, the encounter was not considered for CQMs.
10461	Unexpected Behavior	CCDA	Corrected how the Current Meds Documented checkbox start time is determined and aggregated to the numerator for the measure.
10348	Unexpected Behavior	CPOE	Patient Health Insurance data is now included on the order.
10372	Unexpected Behavior	CQM	Corrected CQM/CDS PowerTab when more than 25 providers were selected.
10410	Unexpected Behavior	CQM	Improved quality of CQM cloud platform to mitigate risk of error stacks being thrown when a CQM report is generated from a client environment.
10466	Unexpected Behavior	CQM	Corrected error in a doctor's CQM report or error stack when entering the CQM/CDS PowerTab if a provider's Shared ID had a trailing space.
10378	Unexpected Behavior	DDE	File generation for various exported file types (i.e., Communicator or CQM) experienced long delays.
10503	Unexpected Behavior	Forms	Improved performance and eliminated time-outs when the v9 Family History Interoperability Dashboard App was used.
10356	Unexpected Behavior	Radiology	Radiology orders were not being generated if the patient did not have a referring provider.
10473	Unexpected Behavior	TOC	Corrected patient matching for Direct Messaging when client used legacy Lab Management Module.

## Features and Fixes for SRS EHR v9.2.071

TTID	Type	Functional Area	Summary
-	New Feature	CQM	Officially certified on 39 CQM measures.
-	New Feature	CPOE	Lab Compendium and Categorization implemented for use with CPOE, Labs, and AMC Reports.
10190	Unexpected Behavior	AMC Report	AMC – Clinical Summary numerator did not get incremented when user had portal access.
8530	Unexpected Behavior	Care Team	Care Team – Providers added from Demographics were not saved under Care Team unless the page was refreshed.
9070	Unexpected Behavior	Clinical Summary	Procedures did not display on the Clinical Summary.
10191	Enhanced Feature	Clinical Summary	Increased font size on the printed copy of the Clinical Summary to improve legibility.
10022	Unexpected Behavior	CPOE	Order Tracking – Matched results text on lab results was cut off when using 1280 x 1024 resolution.
9591	Unexpected Behavior	CPOE	Order Tracking – Did not resize correctly when opened in a smaller SRS window and then maximized.
9045	Unexpected Behavior	CPOE	The Historical Orders window did not resize automatically according to available Interoperability Dashboard screen real estate.
9495	Unexpected Behavior	CPOE	Lab Results – Edit Order Details screen did not resize.
8233	Unexpected Behavior	CPOE	CPOE – Orders Entry screen did not resize if Appointments panel was expanded to the right.
10057	Unexpected Behavior	CQM Report	CMS 68 – When patient had multiple encounters and fulfilled the numerator for each encounter, the CQM report only counted 1 of the encounters in the numerator.
9459	Unexpected Behavior	Dx App	Error stack occurred when clicking the Diagnosis refresh icon more than once on the Interoperability Dashboard.
10286	Unexpected Behavior	Dx App	Diagnosis note was lost when refreshing a diagnosis.
10338	Unexpected Behavior	Dx App	User could not add ICD9/10 codes that contained a double quote within the diagnosis description. (i.e., 361.32).
9983	Unexpected Behavior	Forms	Forms Chart.xml failed when patient had LMM results with NULL or unexpected data.
10324	Unexpected Behavior	Rx Module	RX Drug DB update failed when SQL replication was configured.
10194	Unexpected Behavior	SOP App	MEDICARE (1) did not appear when searched.
10193	Unexpected Behavior	Vitals App	Vitals – Post-screening plan checkbox did not appear when BMI was outside normal range.

## Features and Fixes for SRS EHR v9.1.576

TTID	Type	Functional Area	Summary
-	New Feature	CPOE	Certification for Quest Laboratories interface via Emdeon
-	Enhanced Feature	Forms	Improved system performance for platforms that use Group Tagging.
10248	Unexpected Behavior	AMC Report	Providers were not receiving credit for VDT Measure A after indicating patient consent for SRS Communicator.
10211	Unexpected Behavior	DB Performance	Large HIPAA Audit trail database tables introduced latency when logging into SRS EHR.
10212	Unexpected Behavior	Dx App	Having 2 Diagnosis Super Favorites with the same ICD9 and ICD10 codes caused an SRS App Error.
9969	Enhanced Feature	Forms	Improved the performance of generating the Chart.xml forms data file.



## Features and Fixes for SRS EHR v9.1.566

TTID	Type	Functional Area	Summary
-	Enhanced Feature	CDS	Cloud-based Clinical Decision Support rules.
-	New Feature	Communicator	SRS Communicator Patient Portal Enrollment Security.
-	New Feature	CPOE	Manual Lab Order Entry.
-	New Feature	CPOE	Lab Results – Matched, Unmatched, Solicited, and Unsolicited.
-	New Feature	CPOE	CPOE Future and Completed.
-	New Feature	CPOE	Historical Lab Order Entry.
-	New Feature	CPOE	Lab Order Tracking.
-	New Feature	CPOE	Emdeon Integration.
-	New Feature	CPOE	Emdeon Retrieval.
-	New Feature	CPOE	Radiology Orders and Results integration with HL7-compliant PACs systems.
-	Enhanced Feature	CQM Report	Cloud-based CQM reporting.
-	New Feature	Interoperability Dashboard	Rx Decline Interoperability Dashboard web application.
-	New Feature	Interoperability Dashboard	Consolidated Allergies Interoperability Dashboard web application.
9878	Unexpected Behavior	Communicator	SRS Communicator messages inbound could not be defaulted to a pool.
10179	Unexpected Behavior	Clinical Summary	Clinical Summary displayed incorrect data due to race condition.
9842	Unexpected Behavior	Dx App	Diagnosis Super Favorites caused the web application to fail on load.
9972	Enhanced Feature	Forms	User profile setting introduced to control the creation of the Chart.xml file, among other forms supporting data files.
9990	Enhanced Feature	Intake Manager	Intake Manager respects the SkipChart.xml user profile setting.
9876	Unexpected Behavior	Intake Manager	Missing health insurance plan type caused Intake Manager printing to fail.
9951	Enhanced Feature	Interoperability Dashboard	Improved the performance of the Interoperability Dashboard loading.