

DEPARTMENT OF HEALTH AND HUMAN SERVICES PUBLIC HEALTH SERVICE FOOD AND DRUG ADMINISTRATION GENERAL INFORMATION FIELD TEST RECORD <i>(Use Form FDA 2782, Field Test Record Continuation if more space is needed)</i>		Print Legibly. Use Black Ball Point Pen. Enter one character per box.	FIELD TEST SERIAL NO. (1-8) <div style="border: 1px solid black; padding: 2px; text-align: center; font-weight: bold;">GI</div> REGIONAL REVIEW (NAME)
--	--	--	--


Card No. (9-10)

01	Facility Identification		
02	Name (11-80)		
03	Street Address (11-80)		
04	City (11-73)	State Code <div style="border: 1px solid black; text-align: center; width: 20px; height: 15px; margin: 0 auto;">7475</div>	Zip Code <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">76808184</div>
05	Room Number (11-37)	Person Interviewed (38-56)	Telephone Number <div style="border: 1px solid black; text-align: center; width: 150px; height: 15px; margin: 0 auto;">57596066</div>
Surveyor Information			
06	Name (Last, First, Middle Initial) (11-67)	Accomp. Agency <div style="border: 1px solid black; text-align: center; width: 40px; height: 15px; margin: 0 auto;">6869</div>	Accomp. District <div style="border: 1px solid black; text-align: center; width: 40px; height: 15px; margin: 0 auto;">7072</div>
	Signature	FDA Region <div style="border: 1px solid black; text-align: center; width: 40px; height: 15px; margin: 0 auto;">7374</div>	Date (Mo/Day/Yr) <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">7582</div>
Survey Information			
	Purpose of Survey I -- Initial Survey C -- Compliance F/U H -- HIA R -- Resurvey A -- Audit		
	Assembler Report Nos. (Letter included) <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">121819252632</div>		
	Installation Date Mo/Day/Yr Previous FTRs (5 or 6 digit numbers) <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">334041464752</div>		
System Information			
	Certification Status C -- System Fully Certified V -- Fully Certified With Variance M -- Mixed Certified/Non-Certified N -- Fully Non-Certified		
06	Control Manufacturer	Control Serial No. <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">5263</div>	Date of Mfr. (Mo/Yr) <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">MMYYYY</div>
	Control Model No.	Mfr Code <div style="border: 1px solid black; text-align: center; width: 40px; height: 15px; margin: 0 auto;">6871</div>	Unique ID <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">7277</div>
System Maintenance			
	Is a maintenance schedule designed for compliance being followed? Is the maintenance schedule available for review?		Y -- Yes N -- No X -- Unknown Who does the compliance maintenance? M -- Mfr Rep N -- Not Done P -- Private X -- Unknown I -- In House
06	Instrumentation MDH Serial No. Light Meter Serial No.		
	<div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">11141522</div>		
Number of Each Form Attached		Assembler Data	
06	AR DR CF <div style="border: 1px solid black; text-align: center; width: 20px; height: 15px; margin: 0 auto;">192021</div>		Company
	AF MA HN <div style="border: 1px solid black; text-align: center; width: 20px; height: 15px; margin: 0 auto;">222324</div>		Street Address
	UF KV CT <div style="border: 1px solid black; text-align: center; width: 20px; height: 15px; margin: 0 auto;">252627</div>		City State Zip Code
	MR VC <div style="border: 1px solid black; text-align: center; width: 20px; height: 15px; margin: 0 auto;">282930</div>		Home District FEI Number <div style="border: 1px solid black; text-align: center; width: 100px; height: 15px; margin: 0 auto;">31333443</div>


FORM FDA 3071 (4/01)
EF



spacer assembly
(out of beam)



spacer assembly
(in beam)



beam defining
assembly




6-cm³ ion chamber
(side view)



Al filters



focal spot assembly




slide assembly
(grid down)
with
plastic cassette




test stand
(BENT)



cardboard support
(BENT)



TV view
(slide assembly grid)



BRH test stand



image intensifier
and
spot-film device



0.1 inch Cu



MDH electrometer
(front view)
and
6-cm³ ion chamber



1/8 inch PB

