

Regal-D-Blend®



Mold Chart

for
Anterior and Posterior
Hardened Plastic Teeth

ORTHOCHROMATIC
MULTI-LAYERED
ANTERIOR &
INTERCUSPAL-OCCLUSAL DESIGNED
JUSTI BLEND® POSTERIOR
HARDENED PLASTIC TEETH
SIX SHADE SYSTEM

Highly Cross Linked/Intersectional Polymer Network • Wear Resistant
Vacuum Processed • Fluorescent/Photogenic
No diatoric retention necessary.

100% bondability with methyl methacrylate denture base material.

Complies with ANSI/ADA Specification 15, ISO 3336, BS 3990, NF S91-213, DIN 13914





REGAL-D-BLEND MULTI-BLENDED ANTERIORS 27 UPPERS & 8 LOWERES

A4			4
	10.0 8.3		

D6			6
	12.0 9.0		

A7			7
	10.0 8.8		

E1			1
	9.3 7.3		

B3			3
	9.5 8.0		

E2			2
	10.0 7.8		

B8			8
	11.3 8.8		

E3			3
	10.5 8.3		

C1			1
	10.0 7.5		

E6			6
	10.0 8.5		

C2			2
	10.0 8.0		

E8			8
	11.0 8.8		

C5			5
	11.8 8.8		

F1			1
	8.5 7.8		

C6			6
	10.8 8.5		

F4			4
	9.0 8.5		

D1			1
	10.0 7.5		

G5			5
	10.5 8.8		

D2			2
	10.0 8.0		

G7			7
	10.8 8.8		

D3			3
	10.5 8.3		

H2			2
	10.0 8.0		

Photographs of teeth are within plus or minus .3% of actual size.

H4			4
			51.0 44.67 — —
10.5	8.3		

J2			2
			46.5 40.75 — —
10.0	7.8		

J3			3
			48.0 42.68 — —
10.0	8.0		

J5			5
			53.0 47.14 — —
11.5	9.0		

K4			4
			49.5 43.90 — —
10.8	8.0		

I			C1, D1 E1, F1
			34.5 30.29 — —
9.0	4.8		

2			C2, D2 E2, H2 J2
			36.0 29.86 — —
8.8	4.5		

3			B3, D3 E3, J3
			37.0 30.17 — —
8.3	4.8		

4			A4, F4 H4, K4
			38.0 31.07 — —
9.0	4.8		

5			C5, G5 J5
			40.0 32.80 — —
9.8	5.3		

6			C6, D6 E6
			39.5 33.32 — —
8.5	5.3		

7			A7, G7
			41.0 33.38 — —
9.5	5.0		

8			B8, E8
			41.5 35.77 — —
9.3	5.5		

0° POSTERIORES - JUSTI BLEND 3 UPPERS & LOWERES *

29 MF			6.69
			8.15
			28.5 30.5

31 MF			6.60
			7.62
			30.0 32.8

33 MF			7.31
			8.30
			32.2 34.0

10° POSTERIORES - JUSTI BLEND 3 UPPERS & LOWERES *

30 F			8.41
			9.04
			30.0 33.0

32 F			8.72
			9.63
			31.0 33.5

34 F			9.15
			10.17
			33.0 36.5

20° POSTERIORES - JUSTI BLEND 6 UPPERS & LOWERES *

ANATOMICAL MEDIUM - 33°
JUSTI BLEND 3 UPPERS & LOWERES *



ANATOMICAL LONG - 33°
3 UPPERS & LOWERES



ANTERIOR

A4	Mold Number		4	Suggested Articulation to Lower Anteriors
	Length of Central without Collar			Width of Lower Anteriors on Curve
	Width of Central			Width of Lower Anteriors on Flat
10.0	8.3	10.0 T	49.0 43.78	

POSTERIOR

30 F	Mold Number		8.41	Width 1st Upper Molar			
	Total Width Upper (D)					9.04	Width 1st Lower Molar
	Total Width Lower (E)						

Mold Classification Symbols:

			Combined Symbol
SQUARE	TAPERING	OVOID	

LEGEND:

- D - Straight line measurement of 4 upper posteriors.
- E - Straight line measurement of 4 lower posteriors.

Tooth measurements are within plus or minus 5% of stated value.



ORDER FORM

REGAL-D-BLEND® PRODUCTS:



Order No.: _____ Date: _____
 Bill To: _____
 Address: _____
 City: _____ State: _____ Zip: _____
 Phone # _____ Fax # _____
 Ship To: _____
 Address: _____
 City: _____ State: _____ Zip: _____

Item No.	Description
079-50-0901	Regal-D-Blend Anterior Teeth (1 x 6)
084-50-0903	Regal-D-Blend Maxillary Centrals Shade Guide in 6 Shades
084-50-0904	Regal-D-Blend Anterior Working Mold Guide Contains one set of each mold - 35 Anteriors, all one shade
079-58-0003	Justi Blend Posterior Teeth (1 x 8)
084-59-0011	Justi Blend Posterior Working Mold Guide Contains one set of each mold - 36 posteriors, all one shade

1 x 6 27 UPPER ANTERIOR REGAL-D-BLEND® MOLDS Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
Molds						
A4						
A7						
B3						
B8						
C1						
C2						
C5						
C6						
D1						
D2						
D3						
D6						
E1						
E2						
E3						
E6						
E8						
F1						
F4						
G5						
G7						
H2						
H4						
J2						
J3						
J5						
K4						
Total						

1 x 6 8 LOWER ANTERIOR REGAL-D-BLEND® MOLDS Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
Molds						
1						
2						
3						
4						
5						
6						
7						
8						
Total						

Anterior Cards (1 x 6) - 14 cards per box.

Total Upper Anteriors 1 x 6	
Total Lower Anteriors 1 x 6	
Total Anteriors 1 x 6	

Posterior Cards (1 x 8) - 20 cards per box.

Total Upper Posteriors 1 x 8	
Total Lower Posteriors 1 x 8	
Total Posteriors 1 x 8	

Item #079-50-0901

Item #079-58-0003

Remarks _____

Signed _____

1 x 8 18 UPPER POSTERIOR JUSTI BLEND® MOLDS Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
Molds						
0°						
29MF						
31MF						
33MF						
10°						
30F						
32F						
34F						
20°						
29L						
31L						
33L						
29M						
31M						
33M						
ANATOMICAL MEDIUM - 33°						
30M						
32M						
34M						
ANATOMICAL LONG - 33°						
30L						
32L						
34L						
Total						

1 x 8 18 LOWER POSTERIOR JUSTI BLEND® MOLDS Boxes

Shades	AX/200	AA/202	BB/204	CC/206	DD/208	EE/210
Molds						
0°						
29MF						
31MF						
33MF						
10°						
30F						
32F						
34F						
20°						
29L						
31L						
33L						
29M						
31M						
33M						
ANATOMICAL MEDIUM - 33°						
30M						
32M						
34M						
ANATOMICAL LONG - 33°						
30L						
32L						
34L						
Total						

NOTE: Posteriors are also available in shades 212, 214 and 216 by special request.

This is a sample of our tooth order form. You may use a photo copy of this form for your first order.

See page 6 for details on Justi Blend posterior teeth.

Regal-D-Blend® Articulation Chart

Anteriors		Posteriors			
Upper Molds	Articulates Lower Anteriors	Zero Degree MF	10 Degree F	20 Degree	Anatomical 33 Degree Short, Med., Long
A4	4	31	32	29L	30
A7	7	31	-	33M	32
B3	3	31	32	31M	30
B8	8	31	-	33M	32
C1	1	29	30	29L	30
C2	2	29	30	29L	30
C5	5	31	32	31L	32
C6	6	31	32	31L	32
D1	1	29	30	29M	30
D2	2	31	32	31L	30
D3	3	31	32	29L	30
D6	6	33	34	33L	34
E1	1	29	30	29L	30
E2	2	31	30	29L	30
E3	3	31	32	31L	30
E6	6	31	32	31L	32
E8	8	33	34	33L	32
F1	1	29	30	29L	30
F4	4	31	-	31M	30
G5	5	31	32	31L	32
G7	7	31	32	31L	32
H2	2	31	30	29L	30
H4	4	31	32	31L	30
J2	2	29	32	29L	30
J3	3	29	32	31L	30
J5	5	31	34	33L	32
K4	4	31	-	29L	30



Item# 084-50-0904
Regal-D-Blend Anterior Working Mold Guide.
One set of each mold - 35 anteriors, all one shade.

Item# 084-50-0903
Regal-D-Blend Maxillary Centrals Shade Guide in 6 shades.

* Posterior Teeth

Regal-D-Blend® Anterior Hardened Plastic Teeth have a perfect mate for the posterior arrangement with Justi Blend® Posterior Teeth. The posteriors pictured in this mold chart are the Justi Blend® Posteriors. Justi Blend® Posteriors, like Regal-D-Blend® Anteriors, meet and exceed all ADA Accepted Specifications.

Anterior Mold Comparison With Dentsply's Bioblend®, Trublend®, Portrait®

Regal-D-Blend®	Bioblend® IPN® Bioblend® Trublend®, Portrait®	Regal-D-Blend®	Bioblend® IPN® Bioblend® Trublend®, Portrait®
Uppers		Uppers	
A4	12E	F4	43F
A7	12G	G5	45F
B3	32E	G7	45H
B8	32G	H2	55D
C1	21C	H4	55F
C2	21D	J2	62D
C5	21E	J3	55D
C6	21X	J5	62G
D1	22C	K4	75E
D2	13E	Lowers	
D3	22E	1	C
D6	22H	2	F
E1	42C	3	H
E2	42D	4	N
E3	42F	5	O
E6	42G	6	P
E8	42H	7	R
F1	43D	8	S

Shade Comparison with Major Brands

Color	Regal-D-Blend®	Justi Blend®	Bio-blend®	Bio-form®	Tru-Blend™	Vita®	Por-trait™	Ivoclar Vivadent
White/Gray	AX	200	100	B59	T2	A-1	P-1	-
Yellow/Pink	AA	202	102	B53	T3	A2-B2	P-2	1A
Yellow/Pink	BB	204	104	B92	T4	A2	P-2	1A, 2A
Gray/Pink	CC	206	106	B69	T5	C2-D3	P69	-
Yellow/Brown	DD	208	108	B62	T12	A3-D3	P-3	1C
Gray	EE	210	110	B66/93	T15	A35	P33/P66	1C

This shade comparison chart is to be used in selecting approximate color matches. Slight variations occur due to the difference in the texture and character of the respective materials. This chart should be used as an aid to the skill and preferences of the dentist or technician who must make the final shade choice.

Regal-D-Blend® Anteriors match the central, lateral and cuspid, right & left respectively. Each central, lateral and cuspid are the same shade in three different tones thereby creating a natural blended effect.

Justi Blend® Posteriors match the right or left cuspid body of Regal-D-Blend in tone and shade.

Regal-D-Blend® teeth meet and exceed ADA Specification 15 for Certification.

Distributed by:

Instructions for use of synthetic resin teeth

Note: If you would like to receive written instructions for obtaining optimum results with our teeth please contact:
Justi Products, A Division of American Tooth Industries,
1200 Stellar Drive, Oxnard, CA 93033-2404
Toll Free (800) 235-4639, California (805) 487-9868
Fax (805) 483-8482 • Website: www.americantooth.com

MANUFACTURED BY:



Justi Products

A Brand Division of American Tooth Industries

1200 Stellar Drive, Oxnard, CA 93033-2404

(805) 487-9868 • FAX (805) 483-8482

E-mail: info@americantooth.com • Website: www.americantooth.com