

Luga	Fitter Last Name:			Account Number: 4057807 Patient First Name:
MEDICAL INC.				Fitter First Name:
				(example PT/OT/PTA)
Measure in centin	<u>ieters</u>	_		STOM CIRCULAR KNIT ARM SLEEVE OBS
Left 🗌	Right [
Length Measurem	ent from a	axilla to w	<u>rist (c-g)</u> =	
LONG LENGTH IS >	43 CMS	REGULA	R LENGTH I	S < 43 CMS
				the inside of the arm to the inseam, Axilla , of vith arm down by the side.
Measurement Point	<u>:s</u>			
				n at the axilla. For ease in determining this point, the arm will be even with the top of the book.
	distally or	n the arm v	vhere the w	ession and therefore the most important point. rist meets the hand. One should find the smallest re too tight.
Circumference Wi	<u>ist</u>	=	with s	Silicone with no Silicone
Circumference 'M	<u>id-Lower</u> '	=		
Circumference 'M	id-Upper'	=		
Circumference Ax	<u>killa</u>	=		
Material	Class I 15-20 mmHg	Class II 20-30 mmHg	Class II 30-40 mmHg	Colors Available
Bella Strong				Natural
Bella Lite Arm Sleeve				Beige Brown Black * Silicone only available in Beige * Brown and Black are only available in 20-30 mmHg
Sizing Chart - Bella Lite Armsleev	re			Sizing Chart - Bella Strong Armsleeve
				Size Wrist Mid forearm Axilla (C) (D) (G)
				1 14 18 cm 17 - 21 cm 22 - 31 cm
				2 14-18cm 17-21cm 31-41cm 3 14-18cm 21-25cm 25-34cm
Size Wrist	Mid-Low	or	Axilla	4 16-19cm 21-25cm 28-38cm 5 17-20cm 21-25cm 34-43cm
Size Wilst	W.C.			6 14-18cm 25-29cm 29-38cm
c 5%"-75%" (14cm - 20cm)	-7"- 9½" (18cm -	23cm) 8½"	15½" (22cm - 40	7 17-20 cm 25-29 cm 33-43 cm 8 18-21 cm 25-29 cm 38-48 cm
3	/ - 978 (16CH) - /4"- 111/4" (23cm		- 17 ³ / ₄ " (27cm - 40	9 16-19cm 29-33cm 32-42cm
L 6 ⁷ /8" - 9 ¹ /4" (18cm - 24cm) 11 ¹			- 19 ⁷ /8" (33cm - 5	
XL				All sizes are available in 2 lengths:
				Regular (up to 43cm) Long (greater than 43cm)