Order Form

Datamatrix Coding Service for Cryo.s™

Cryo.s™ with Datamatrix Code

- S ECC 200 Datamatrix
- Laser written Datamatrix code with a strong resistance to chemical and mechanical damage
- Minimum order quantity: 1 case (equivalent to 500 1 ml or 2 ml Cryo.s™ or 300 4 ml Cryo.s™)
- 100% verification of the readability of each Datamatrix
- Encrypted information may represent sequential sample IDs, with the option of prefixed alphanumeric characters (e.g. ABC00001-ABC10000), or consistent alphanumeric characters (e.g. CHO-K1 on each tube)

High data capacity:

Symbol size	Data capacity (No. of characters)		
(Rows x columns)	numeric	alphanumeric	
10 x 10	6	3	
12 x 12	10	6	
14 x 14	16	10	

Customer No.	Phone
Name	Fax
Institution	E-Mail
Street	Order No. Order Date
City / Zip Code	Shipping Address(if different than above)
Country	

1. Please choose the appropriate Cryo.s™:



☐ Cat.-No. 123 2XX Cryo.s TM , 1 ml, with internal thread and starfoot base



☐ Cat.-No. 126 2XX Cryo.s™, 2 ml, with external thread and starfoot base



☐ Cat.-No. 122 2XX Cryo.s[™], 2 ml, with internal thread and starfoot base



☐ Cat.-No. 127 2XX Cryo.s™, 4 ml, with external thread and starfoot base

□ green

2. Please choose the appropriate screw cap colour:

natura
red

blue
yellow

3. Please choose the desired symbol size:

11	\cap	V	1	\cap

□ 12 x 12

□ 14 x 14

4. Please choose between a sequential series of sample IDs or consistent characters (fixed code):

☐ Sequential series of sample IDs

Sta	rt	of	sequence
_		_	

(e.g. ABC0000001)

End of sequence

(e.g. ABC0010000)

Number of characters

(e.g. 10)

☐ Fixed code

(e.g. CHO-K1 Clone 123)

5. Please specify the required amount of Cryo.s™ with Datamatrix:

	-	_	-	_	
(a a	1	Λ	Λ	\cap	۱۱

Released by (signature)

Place, date

_	_	_	
ī			