

Special GC Offer



Buy a GC package at list price and receive a free Auto Injector or free Auto Injector/Auto Sampler

Combining ease of use with extraordinary speed and precision, the **GC-2010** offers enhanced productivity and performance flexibility for a wide range of applications. It features up to 4 high-sensitivity detectors, a large oven and superior Advanced Flow Control (AFC) design, enabling excellent retention time reproducibility and optimization of high-pressure injections.

Ideal for routine analysis, the **GC-2014** utilizes the highend technology of the GC-2010 to cover a wide range of applications from capillary to packed column analysis.



Increase laboratory productivity by automating your analyses with Shimadzu's AOC-20 sample injection system. Proven injection reproducibility is combined with easy operation and high throughput of up to 150 samples.

- AOC-20i Auto Injector Features enhanced needle guide and injection sequence management, as well as a variety of injection modes
- A0C-20s Auto Sampler Improves productivity with self-correcting software, automatic identification of ex changed racks, and the ability to program for subsequent use in other applications



- 1) Free AOC 20i with single injector/detector model
- 2) Free AOC 20i and S with multiple injectors and detectors

Please reference the following part numbers:

- 220-94708-01: AOC-20i for GC-2010
- 220-94708-02: AOC-20i for GC-2014
- 220-94708-03: AOC-20s

Valid in SSI international regions only. Not valid with any other offer. Listed part #s must be used. Expires 9/30/2009.